

## Laboratory automation systems and workcells

<p><b>Part 1 of 14</b></p> <p><i>See accompanying article on page 12</i>  <i>See captodayonline.com/productguides for an interactive version of guide</i></p>	<p><b>Abbott Diagnostics</b>          Deborah Anderson deborah.anderson@abbott.com          100 Abbott Park Road, Abbott Park, IL 60064          847-936-6353 www.abbottdiagnostics.com</p>	<p><b>Aim Lab (formerly Ai Scientific)</b>          Ralph Donaldson aimlab@aimlab.com          10-22 Hornibrook Esplanade, Clontarf, QL, Australia 4035          +61 7 3105 5005 www.aimlab.com</p>
<p><b>Name of system/First year installed/No. of 2011 contracts signed</b>  <b>No. of live sites installed in N. America/Europe/Asia-Australia</b></p>	<p><b>Accelerator APS/2005/—</b>  <b>20+/-/—</b></p>	<p><b>PathFinder 350S/2008/27</b>  <b>~14 (North and South America)/~10/8</b></p>
<p><b>Automation products that are available</b></p> <ul style="list-style-type: none"> <li>• Pre-analytical processor/Total laboratory automation</li> <li>• Automated functions: Accessioning/Track load/Centrifugation/Decapping</li> <li>• Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing</li> <li>• Automated functions: Storage-retrieval/Intelligent sample routing</li> <li>• SW: Dedicated Process Control/Middleware control using LIS/Architecture</li> <li>• Company has dedicated automation support team/Remote system monitoring</li> </ul>	<p>yes/yes          yes/yes/yes/yes          yes/no/no/yes          yes/yes          yes/yes/open          yes/yes</p>	<p>yes/no          yes/yes/no/yes          yes/yes/yes/yes          no/yes          no/no/open          yes/yes</p>
<p><b>Software features/functionality</b></p> <ul style="list-style-type: none"> <li>• Patient demographics and insurance data/Rules-based architecture</li> <li>• Supports data retrieval/Internet connectivity</li> <li>• Online real-time help system/QC/Stats and management reports</li> <li>• Evaluates validity and releasability of results from automated analyzers</li> <li>• Specimen tracking/Priority processing/Random-access spec. movement</li> <li>• Supports accession number redundancy (duplicate specimen ID)</li> <li>• Supports specimen carrier and level identification</li> <li>• Unique bar-code number per container required</li> <li>• Specimen routing/Multistop routing (one tube to multiple workstations)</li> <li>• Specimen scheduling/Instrument scheduling</li> <li>• Routes test to workstation/Automatic reflex, repeat, dilutions</li> <li>• Supports multiple HW configuration/Supports other proprietary transport. HW</li> <li>• Sample storage and retrieval SW/Supports approved CLSI standards</li> </ul>	<p>automation SW feature/automation SW feature          automation SW feature/automation SW feature          automation SW feature/automation SW feature/automation SW feature          automation SW feature          automation SW feature/automation SW feature/automation SW feature          automation SW feature          automation SW feature          automation SW feature/automation SW feature          automation SW feature/automation SW feature          automation SW feature/automation SW feature          automation SW feature/automation SW feature          automation SW feature/automation SW feature          automation SW feature/automation SW feature          automation SW feature/automation SW feature</p>	<p>—          —/automation SW feature          —/—/automation SW feature          LIS feature          automation SW feature/automation SW feature/—          automation SW feature          automation SW feature          —          automation SW feature/automation SW feature          —          automation SW feature/automation SW feature          automation SW feature/—          LIS feature/automation SW feature</p>
<p><b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b></p>	<p>Cerner Classic, Cerner Millennium, Cortex, Delphic, Dianoema, GE Ultra, GLMIS by MIPS, Lab Track, Medisolution by Technidata, Meditech 5.4, Misys, Misys CPR (Cloverleaf Engine), Misys Smart, ModulabGold (Izasa), OSM, Roche Omega, SCS, Siemens, Soft/HL7, ASTM</p>	<p>Instrument Manager, Ultra, Apollo, Kestral, others/ASTM, CLSI-LIS2A</p>
<p><b>Transportation systems available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Containers device accommodates/Average throughput in cm per second</li> <li>• Supports automatic rerouting for reflex-repeat-dilutions</li> <li>• Modular HW/Installed options/Device can operate in track and manual mode</li> <li>• Required utilities/Required maintenance</li> <li>• Carrier type/Scalable system</li> </ul>	<p>yes          APS track section/40.2 × variable × 17.0 inches/yes          16, 13 × 100; 16, 13 × 75, others, multiple types simultaneously/13          yes          yes/floor mounted/yes          compressed air, electricity, water/—          single specimen container per carrier/yes</p>	<p>no          —          —          —          —          —</p>
<p><b>Automated centrifugation available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Maximum throughput/Containers device accommodates</li> <li>• Can identify tube types for custom programmed rate and spin times per run</li> <li>• More than one centrifuge can be connected to track system</li> <li>• For multi-unit centrifuge, each centrifuge operates independently for rate and time</li> <li>• Maintenance required</li> </ul> <p><b>Automated input/accessioning available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Dedicated lanes for stat samples</li> <li>• Maximum No. of samples that can be loaded/Maintenance required</li> </ul> <p><b>Automated decapping available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> <li>• Removes multiple size tube caps per run/Removes screw type sample caps</li> </ul> <p><b>Automated sorting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Software can sort by</li> </ul> <p><b>Specimen integrity monitor available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> </ul> <p><b>Automated aliquotting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates</li> <li>• Inspects samples for bar code/Detects and reports clots in specimen</li> <li>• Detects and reports quantity not sufficient specimens/Maintenance required</li> </ul>	<p>yes          Hettich/58.5 × 32 × 42 inches/yes          up to 320/16, 13×100; 16, 13×75, others, multiple types simultaneously          no          yes          no          weekly, monthly          yes          input-output module/54.3 × 77.6 × 39.6 inches/yes/up to 600          16, 13 × 100; 16, 13 × 75, others, multiple types simultaneously/yes          744/weekly, monthly          yes          decapper module/46.7 × 34.7 × 17 inches/yes/up to 600          16, 13 × 100; 16, 13 × 75, others, multiple types simultaneously/daily, weekly          yes/yes          yes          input output module/54.3 × 77.6 × 39.6 inches/yes/up to 600          16, 13 × 100; 16, 13 × 75, others, mult. types simult./specimen, method, output          no          —          —          no          —          —          —          —          —</p>	<p>no          —          —          —          —          —          yes          PathFinder 350S/98 × 40 × 52 cm (39 × 16 × 21 inches)/yes/350 tubes per hour          16, 13 × 100; 16, 13 × 75/yes          flexible/annually          no          —          —          —          yes          PathFinder 350S/98 × 40 × 52 cm (39 × 16 × 21 inches)/yes/350 tubes per hour          16, 13 × 100; 16, 13 × 75/specimen, method, output          no          —          —          no          —          —          —</p>
<p><b>Instrument (analyzer) interfaces</b></p> <ul style="list-style-type: none"> <li>• Rules-based instrument interface control subsystem</li> <li>• Process control of instrument via control subsystem</li> </ul> <p><b>Physical/hardware (instrument/specimen) interface</b></p> <ul style="list-style-type: none"> <li>• Hematology/Chemistry/Coagulation</li> <li>• Immunoassay/Urinalysis</li> </ul>	<p>yes          yes          no/point-of-reference sampling/no          point-of-reference sampling/no</p>	<p>no          no          —          —</p>
<p><b>Instruments to which your system or product is interfaced</b>  <b>Other robotic products/components to which system or product is linked</b></p>	<p>Architect c8000, c16000, i2000SR, Ortho Fusion 5.1, Diasorin Liaison (ex-US only)          —</p>	<p>—          —</p>
<p><b>Automated recapper or sealer available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Recaps-seals multiple size tubes simult./Containers device accommodates</li> <li>• Maintenance required</li> </ul>	<p>sealer          resealer module/49.2 × 44.9 × 17 inches/yes/up to 600          yes/16, 13 × 100; 16, 13 × 75, others          monthly</p>	<p>no          —          —          —</p>
<p><b>Automated storage and retrieval available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Connects to the track</li> <li>• Room temperature/Minimum and maximum number of tubes stored per module</li> <li>• Multiple size tubes can be stored in the same module/Maintenance required</li> <li>• Refrigerated storage and retrieval capability</li> <li>• Longitudinal upgrade pathway or plan to protect users' investments</li> <li>• Average time to install/Who provides service, support/Hours support is available</li> <li>• On-site biomedical engineer required/User group meets regularly</li> </ul>	<p>yes          tube storage module/95 × 89.2 × 70 inches/yes/up to 600          16, 13 × 100; 16, 13 × 75, others, multiple types simultaneously/yes          no/0 and 15,360          yes/daily, monthly          yes          modular open architecture          depends on configuration/Abbott Diagnostics/business and extended hours          yes/yes</p>	<p>yes          PathFinder 350S/52 × 98 × 40 cm/yes/350+          16, 13 × 100; 16, 13 × 75/yes          yes/250          yes/weekly, six months          no          ability to network multiple instruments          1 day/distributor/—          no/no</p>
<p><b>List price</b></p> <p><b>Individual list prices for components</b></p> <ul style="list-style-type: none"> <li>• Process control SW/Transportation systems/Auto. centrifugation</li> <li>• Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval</li> <li>• Specimen integrity monitor/Automated aliquot</li> <li>• Instrument (analyzer) interfaces/Automated recap</li> </ul>	<p>varies          —          —          —          —</p>	<p>\$A62,000          included/—/—          included/—/included/—          —          —</p>
<p><b>Distinguishing features (supplied by company)</b>  <i>* For basic bulding block unit</i>  <i>** Average throughput in specimen containers per hour per device</i>  <i>Note: a dash in lieu of an answer means company did not answer question or question is not applicable</i></p>	<p>flexibility: configurable, component-based design; functionality:          refrigerated online storage and multiple tube types simultaneously, RFID,          point-in-space sampling; support: Lean Six Sigma Black Belt consultants;          Class 2 laser; operations manual on Web site</p>	<p>benchtop sorting at an affordable price; ability to change deck layout in one          minute; flexible input and output positions</p>

## Laboratory automation systems and workcells

<b>Part 2 of 14</b>	<b>Aim Lab (formerly Ai Scientific)</b> Ralph Donaldson sales@aimlab.com 10-22 Hornibrook Esplanade, Clontarf, QL, Australia 4035 +61 7 3105 5005 www.aimlab.com	<b>Beckman Coulter</b> Mike Hoang mbhoang@beckman.com 200 S. Kraemer Boulevard, Brea, CA 92822 714-961-6385 www.beckmancoulter.com
See captodayonline.com/productguides for an interactive version of guide		
<b>Name of system/First year installed/No. of 2011 contracts signed No. of live sites installed in N. America/Europe/Asia-Australia</b>	<b>PathFinder 900/2008/7 0/12/7</b>	<b>AutoMate 800/2006/21 30/100/15</b>
<b>Automation products that are available</b> • Pre-analytical processor/Total laboratory automation • Automated functions: Accessioning/Track load/Centrifugation/Decapping • Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing • Automated functions: Storage-retrieval/Intelligent sample routing • SW: Dedicated Process Control/Middleware control using LIS/Architecture • Company has dedicated automation support team/Remote system monitoring	yes/no yes/yes/no/yes yes/yes/yes/yes no/yes no/no/open yes/yes	yes/no yes/no/yes/yes yes/yes/yes/no yes/yes yes/no/open yes/—
<b>Software features/functionality</b> • Patient demographics and insurance data/Rules-based architecture • Supports data retrieval/Internet connectivity • Online real-time help system/QC/Stats and management reports • Evaluates validity and releasability of results from automated analyzers • Specimen tracking/Priority processing/Random-access spec. movement • Supports accession number redundancy (duplicate specimen ID) • Supports specimen carrier and level identification • Unique bar-code number per container required • Specimen routing/Multistop routing (one tube to multiple workstations) • Specimen scheduling/Instrument scheduling • Routes test to workstation/Automatic reflex, repeat, dilutions • Supports multiple HW configuration/Supports other proprietary transport. HW • Sample storage and retrieval SW/Supports approved CLSI standards	—/LIS feature automation SW feature/automation SW feature automation SW feature/—/automation SW feature LIS feature automation SW feature/automation SW feature/automation SW feature automation SW feature automation SW feature — automation SW feature/automation SW feature LIS feature/LIS feature LIS feature/LIS feature automation SW feature/— LIS feature/automation SW feature	LIS feature/automation SW feature LIS feature/— automation SW feature/LIS feature/automation SW feature LIS feature automation SW feature/automation SW feature/automation SW feature automation SW feature automation SW feature automation SW feature automation SW feature/automation SW feature automation SW feature/— automation SW feature/— automation SW feature/— automation SW feature/automation SW feature
<b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b>	Instrument Manager, Ultra, Kestral, Apollo, others//ASTM, CLSI-LIS2A	SCC, Siemens, Philips/ASTM, Power Processor
<b>Transportation systems available</b> • Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Containers device accommodates/Average throughput in cm per second • Supports automatic rerouting for reflex-repeat-dilutions • Modular HW/Installed options/Device can operate in track and manual mode • Required utilities/Required maintenance • Carrier type/Scalable system	no — — — — — —	no — — — — — —
<b>Automated centrifugation available</b> • Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Maximum throughput/Containers device accommodates • Can identify tube types for custom programmed rate and spin times per run • More than one centrifuge can be connected to track system • For multi-unit centrifuge, each centrifuge operates independently for rate and time • Maintenance required	no — — — — — —	yes AutoMate 800/—/yes 300/16, 13 × 100; 16, 13 × 75, Sarstedt, Greiner, BD pediatric tubes no no no daily yes
<b>Automated input/accessioning available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Dedicated lanes for stat samples • Maximum No. of samples that can be loaded/Maintenance required	yes PathFinder 900 module/1.7 × 2.5 × 1.4 m/yes/350 tubes per hour 16, 13 × 100; 16, 13 × 75/yes up to 250/quarterly	yes AutoMate 800/—/yes/420 16, 13 × 100; 16, 13 × 75, Sarstedt, Greiner, BD pediatric tubes/yes 600/daily, monthly
<b>Automated decapping available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required • Removes multiple size tube caps per run/Removes screw type sample caps	yes PathFinder 900 module/1.7 × 2.5 × 1.4 m/yes/>500 tubes per hour 16, 13 × 100; 16, 13 × 75/weekly, monthly, annually yes/yes	yes AutoMate 800/—/yes/420 16, 13 × 100; 16, 13 × 75, Sarstedt, Greiner, BD pediatric/daily, monthly yes/yes
<b>Automated sorting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Software can sort by	yes PathFinder 900 module/1.7 × 2.5 × 1.4 m/yes/350 tubes per hour 16, 13 × 100; 16, 13 × 75/specimen, method, output	yes AutoMate 800/—/yes/420 16, 13 × 100; 16, 13 × 75, Sarstedt, Greiner, BD pediatric/method, output
<b>Specimen integrity monitor available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required	no — —	no — —
<b>Automated aliquotting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates • Inspects samples for bar code/Detects and reports clots in specimen • Detects and reports quantity not sufficient specimens/Maintenance required	yes PathFinder 900 module/1.7 × 2.5 × 1.4 m/yes/>500 tubes per hour 16, 13 × 100; 16, 13 × 75 yes/yes yes/weekly, monthly, annually	yes AutoMate 800/—/yes/420 16, 13 × 100; 16, 13 × 75, Sarstedt yes/yes yes/daily, monthly
<b>Instrument (analyzer) interfaces</b> • Rules-based instrument interface control subsystem • Process control of instrument via control subsystem Physical/hardware (instrument/specimen) interface • Hematology/Chemistry/Coagulation • Immunoassay/Urinalysis	no no no/no/no no/no	no no — —
<b>Instruments to which your system or product is interfaced Other robotic products/components to which system or product is linked</b>	none —	— —
<b>Automated recapper or sealer available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Recaps-seals multiple size tubes simult./Containers device accommodates • Maintenance required	yes PathFinder 900 module/1.7 × 2.5 × 1.4 m/yes/>600 tubes per hour no/16, 13 × 100; 16, 13 × 75 daily, monthly, annually	no — — —
<b>Automated storage and retrieval available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Connects to the track • Room temperature/Minimum and maximum number of tubes stored per module • Multiple size tubes can be stored in the same module/Maintenance required • Refrigerated storage and retrieval capability Longitudinal upgrade pathway or plan to protect users' investments	yes PathFinder 900/1.7 × 2.5 × 1.4 m/yes/900+ 16, 13 × 100; 16, 13 × 75/yes no/1,000 in standard format yes/weekly, monthly, annually no ability to network multiple PathFinders	yes AutoMate 800/—/yes/420 16, 13 × 100; 16, 13 × 75, Sarstedt, Greiner, BD pediatric tubes/no yes/1 and 400 yes/daily, monthly no —
<b>Average time to install/Who provides service, support/Hours support is available On-site biomedical engineer required/User group meets regularly</b>	3 weeks/GST and distributor/24-7 no/no	7 days/Beckman Coulter/24-7 no/no
<b>List price</b>	\$A420,000 (fully optioned system)	—
<b>Individual list prices for components</b> • Process control SW/Transportation systems/Auto. centrifugation • Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval • Specimen integrity monitor/Automated aliquot • Instrument (analyzer) interfaces/Automated recap	included/included/— included/included/included/— —/included —/included	— — — —
<b>Distinguishing features (supplied by company)</b>	ability to change deck layout in five minutes; ability to recap primary tubes with original cap to preserve tube-type identity; dual track for parallel processing of samples leading to high overall throughput (that is, output tubes)	automatic rack layout can be reconfigured with another rack style; intelligent aliquotting; sample storage routing by duration and temperature
* For basic building block unit ** Average throughput in specimen containers per hour per device Note: a dash in lieu of an answer means company did not answer question or question is not applicable		

## Laboratory automation systems and workcells

<p><b>Part 3 of 14</b></p> <p>See <a href="http://captodayonline.com/productguides">captodayonline.com/productguides</a> for an interactive version of guide</p>	<p><b>Beckman Coulter</b>                  Christian Frenz cfrenz@beckman.com                  200 S. Kraemer Boulevard, Brea, CA 92822                  714-961-6385 www.beckmancoulter.com</p>	<p><b>Beckman Coulter</b>                  Mike Hoang mbhoang@beckman.com                  200 S. Kraemer Boulevard, Brea, CA 92822                  714-961-6385 www.beckmancoulter.com</p>
<p><b>Name of system/First year installed/No. of 2011 contracts signed</b>  <b>No. of live sites installed in N. America/Europe/Asia-Australia</b></p>	<p><b>AutoMate 2500 Family/2003/85</b>                  77/485/42</p>	<p><b>LH 1500 Hematology Automation Series/2002/6</b>                  100/20/20</p>
<p><b>Automation products that are available</b></p> <ul style="list-style-type: none"> <li>• Pre-analytical processor/Total laboratory automation</li> <li>• Automated functions: Accessioning/Track load/Centrifugation/Decapping</li> <li>• Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing</li> <li>• Automated functions: Storage-retrieval/Intelligent sample routing</li> <li>• SW: Dedicated Process Control/Middleware control using LIS/Architecture</li> <li>• Company has dedicated automation support team/Remote system monitoring</li> </ul>	<p>yes/no                  yes/no/no/yes                  yes/yes/yes/yes                  no/yes                  yes/yes/open                  yes/yes</p>	<p>yes/yes                  yes/yes/no/no                  yes/no/no/no                  yes/yes                  yes/yes/open                  yes/yes</p>
<p><b>Software features/functionality</b></p> <ul style="list-style-type: none"> <li>• Patient demographics and insurance data/Rules-based architecture</li> <li>• Supports data retrieval/Internet connectivity</li> <li>• Online real-time help system/QC/Stats and management reports</li> <li>• Evaluates validity and releasability of results from automated analyzers</li> <li>• Specimen tracking/Priority processing/Random-access spec. movement</li> <li>• Supports accession number redundancy (duplicate specimen ID)</li> <li>• Supports specimen carrier and level identification</li> <li>• Unique bar-code number per container required</li> <li>• Specimen routing/Multistop routing (one tube to multiple workstations)</li> <li>• Specimen scheduling/Instrument scheduling</li> <li>• Routes test to workstation/Automatic reflex, repeat, dilutions</li> <li>• Supports multiple HW configuration/Supports other proprietary transport. HW</li> <li>• Sample storage and retrieval SW/Supports approved CLSI standards</li> </ul>	<p>LIS feature/automation SW feature                  automation SW feature/—                  automation SW feature/—/automation SW feature                  —                  automation SW feature/automation SW feature/automation SW feature                  automation SW feature                  automation SW feature                  —                  automation SW feature/automation SW feature                  automation SW feature/—                  —/—                  automation SW feature/—                  automation SW feature/automation SW feature</p>	<p>—                  automation SW feature/—                  automation SW feature/LIS feature/—                  —                  automation SW feature/automation SW feature/automation SW feature                  —                  —                  automation SW feature                  automation SW feature/ automation SW feature                  automation SW feature/ automation SW feature                  automation SW feature/ automation SW feature                  automation SW feature/—                  automation SW feature/ automation SW feature</p>
<p><b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b></p>	<p>Cerner, Misys, Modulus, Data Innovations, SCC, Atlas, McKesson/HL7, ASTM</p>	<p>Cerner, Sunquest, SCC, Meditech, others/LH 1500</p>
<p><b>Transportation systems available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions* (H x W x D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Containers device accommodates/Average throughput in cm per second</li> <li>• Supports automatic rerouting for reflex-repeat-dilutions</li> <li>• Modular HW/Installed options/Device can operate in track and manual mode</li> <li>• Required utilities/Required maintenance</li> <li>• Carrier type/Scalable system</li> </ul>	<p>no                  —                  —                  —                  —                  —</p>	<p>yes                  —/—/yes                  13 x 75/—                  yes                  yes/floor mounted/yes                  compressed air, electricity/monthly                  single specimen container per carrier/yes</p>
<p><b>Automated centrifugation available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions (H x W x D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Maximum throughput/Containers device accommodates</li> <li>• Can identify tube types for custom programmed rate and spin times per run</li> <li>• More than one centrifuge can be connected to track system</li> <li>• For multi-unit centrifuge, each centrifuge operates independently for rate and time</li> <li>• Maintenance required</li> </ul> <p><b>Automated input/accessioning available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Dedicated lanes for stat samples</li> <li>• Maximum No. of samples that can be loaded/Maintenance required</li> </ul> <p><b>Automated decapping available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> <li>• Removes multiple size tube caps per run/Removes screw type sample caps</li> </ul> <p><b>Automated sorting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Software can sort by</li> </ul> <p><b>Specimen integrity monitor available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> </ul> <p><b>Automated aliquotting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates</li> <li>• Inspects samples for bar code/Detects and reports clots in specimen</li> <li>• Detects and reports quantity not sufficient specimens/Maintenance required</li> </ul>	<p>no                  —                  —                  —                  —                  —                  —                  yes                  AutoMate 2500 Family units/64 x 73 x 53 inches/yes/1,200                  16, 13 x 100; 16, 13 x 75; diameter: 10.5-17.0 mm; length: 70-100 mm/yes                  300, continuously/—                  —                  yes                  AutoMate 2500 Family units/64 x 73 x 53 inches/yes/1,200                  16, 13 x 100; 16, 13 x 75; diameter: 10.5-17.0 mm; length: 70-100 mm/—                  yes/yes                  —                  yes                  AutoMate 2500 Family units/64 x 73 x 53 inches/yes/1,200                  16, 13x100; 16, 13x75; others/specimen, test order, fill level, input position                  no                  —                  —                  yes                  AutoMate 1250, 2550/64 x 101 x 53 inches/yes/600                  16, 13 x 100; 16, 13 x 75, secondary tubes 13 x 75                  yes/yes                  yes/daily</p>	<p>no                  —                  —                  —                  —                  —                  —                  yes                  —                  13 x 75/yes                  200/monthly                  no                  —                  —                  —                  yes                  —/—/yes/425                  13 x 75/method                  no                  —                  —                  no                  —                  —                  —</p>
<p><b>Instrument (analyzer) interfaces</b></p> <ul style="list-style-type: none"> <li>• Rules-based instrument interface control subsystem</li> <li>• Process control of instrument via control subsystem</li> </ul> <p><b>Physical/hardware (instrument/specimen) interface</b></p> <ul style="list-style-type: none"> <li>• Hematology/Chemistry/Coagulation</li> <li>• Immunoassay/Urinalysis</li> </ul>	<p>no                  no                  no/no/no                  no/no</p>	<p>no                  yes                  robotic arm interface/—/—                  —</p>
<p><b>Instruments to which your system or product is interfaced</b>  <b>Other robotic products/components to which system or product is linked</b></p>	<p>—                  —</p>	<p>LH 750, 755 and LH 780, 785                  —</p>
<p><b>Automated recapper or sealer available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput*</li> <li>• Recaps-seals multiple size tubes simult./Containers device accommodates</li> <li>• Maintenance required</li> </ul>	<p>sealer                  all AutoMate 2500 Family units/—/yes/1,200                  yes/16, 13 x 100; 16, 13 x 75                  daily</p>	<p>no                  —                  —                  —</p>
<p><b>Automated storage and retrieval available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput*</li> <li>• Containers device accommodates/Connects to the track</li> <li>• Room temperature/Minimum and maximum number of tubes stored per module</li> <li>• Multiple size tubes can be stored in the same module/Maintenance required</li> <li>• Refrigerated storage and retrieval capability</li> </ul> <p><b>Longitudinal upgrade pathway or plan to protect users' investments</b></p>	<p>—                  —                  —                  —                  —                  longitudinal upgrade pathway</p>	<p>yes                  —/—/yes/340                  13 x 75/yes                  yes/1,000                  no/weekly, monthly                  —                  expandable, as the lab grows</p>
<p><b>Average time to install/Who provides service, support/Hours support is available</b>  <b>On-site biomedical engineer required/User group meets regularly</b></p>	<p>1 week/Beckman Coulter/24-7                  no/no</p>	<p>7-21 days/Beckman Coulter/24-7                  no/yes</p>
<p><b>List price</b></p>	<p>\$290,000-\$460,000</p>	<p>varies</p>
<p><b>Individual list prices for components</b></p> <ul style="list-style-type: none"> <li>• Process control SW/Transportation systems/Auto. centrifugation</li> <li>• Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval</li> <li>• Specimen integrity monitor/Automated aliquot</li> <li>• Instrument (analyzer) interfaces/Automated recap</li> </ul>	<p>—                  —                  —                  —</p>	<p>—                  —                  —                  —</p>
<p><b>Distinguishing features (supplied by company)</b></p> <p>* For basic building block unit                  ** Average throughput in specimen containers per hour per device                  Note: a dash in lieu of an answer means company did not answer question or question is not applicable</p>	<p>high speed, cost efficient way to automate pre and postanalytical steps; improves patient safety and lab efficiency through tube inspection unit to ensure correct label is on the sample and that enough sample volume is available; allows direct sorting to most analyzers' racks, and easy to change configurations</p>	<p>automatic hands-off rerun and reflex test from the stockyard to the analyzers; sorting of pending samples for secondary tests by test; automatically loads analyzers and is expandable</p>

Tabulation does not represent an endorsement by the College of American Pathologists.



## Laboratory automation systems and workcells

# Online Now

**More  
instrument  
and  
software  
system  
product  
guides  
online**

## CAP TODAY INTERACTIVE PRODUCT GUIDES

Coagulation analyzers

Hematology analyzers

Laboratory information systems

Next-generation sequencing

Chemistry analyzers

Blood bank  
information systems

In vitro blood gas analyzers

Automated molecular platforms

Positive patient ID

Automated immunoassay

POC coagulation analyzers

Billing/accounts  
receivable systems

Laboratory-provider  
links software

Bedside glucose  
testing systems

Middleware systems

Anatomic pathology  
computer systems

**GO TO:**

[captodayonline.com/  
productguides](http://captodayonline.com/productguides)

Part 4 of 14

See [captodayonline.com/productguides](http://captodayonline.com/productguides) for an interactive version of guide

Beckman Coulter

Mike Hoang [mbhoang@beckman.com](mailto:mbhoang@beckman.com)  
200 S. Kraemer Boulevard, Brea, CA 92822  
714-961-6385 [www.beckmancoulter.com](http://www.beckmancoulter.com)

Name of system/First year installed/No. of 2011 contracts signed  
No. of live sites installed in N. America/Europe/Asia-Australia

Power Processor/1998/38  
402/126/146

Automation products that are available

- Pre-analytical processor/Total laboratory automation
- Automated functions: Accessioning/Track load/Centrifugation/Decapping
- Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing
- Automated functions: Storage-retrieval/Intelligent sample routing
- SW: Dedicated Process Control/Middleware control using LIS/Architecture
- Company has dedicated automation support team/Remote system monitoring

yes/yes  
yes/yes/yes/yes  
yes/yes/yes/yes  
yes/yes  
yes/yes/open  
yes/yes

Software features/functionality

- Patient demographics and insurance data/Rules-based architecture
- Supports data retrieval/Internet connectivity
- Online real-time help system/QC/Stats and management reports
- Evaluates validity and releasability of results from automated analyzers
- Specimen tracking/Priority processing/Random-access spec. movement
- Supports accession number redundancy (duplicate specimen ID)
- Supports specimen carrier and level identification
- Unique bar-code number per container required
- Specimen routing/Multistop routing (one tube to multiple workstations)
- Specimen scheduling/Instrument scheduling
- Routes test to workstation/Automatic reflex, repeat, dilutions
- Supports multiple HW configuration/Supports other proprietary transport. HW
- Sample storage and retrieval SW/Supports approved CLSI standards

LIS feature/automation SW feature  
automation SW feature/automation SW feature  
automation SW feature/automation SW feature/automation SW feature  
automation SW feature  
automation SW feature/automation SW feature/automation SW feature  
—  
automation SW feature  
automation SW feature  
automation SW feature/automation SW feature  
automation SW feature/automation SW feature  
automation SW feature/automation SW feature  
automation SW feature/—  
automation SW feature/automation SW feature

LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS

SCC, Siemens, Philips, Misys, Cerner, McKesson, GE, Mediatech, PerSe, Molis, MIPS, Vista, Swiss Lab/Power Processor, Direct, HL7

Transportation systems available

- Model/Dimensions\* (H × W × D)/Conforms to CLSI Stand. Auto 1-5
- Containers device accommodates/Average throughput in cm per second
- Supports automatic rerouting for reflex-repeat-dilutions
- Modular HW/Installed options/Device can operate in track and manual mode
- Required utilities/Required maintenance
- Carrier type/Scalable system

yes  
Power Processor II/—/yes  
16, 13 × 100; 16, 13 × 75, Sarstedt/—  
yes  
yes/floor and subfloor mounted/yes  
compressed air, electricity/monthly  
single specimen container per carrier/yes

Automated centrifugation available

- Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5
- Maximum throughput/Containers device accommodates
- Can identify tube types for custom programmed rate and spin times per run
- More than one centrifuge can be connected to track system
- For multi-unit centrifuge, each centrifuge operates independently for rate and time
- Maintenance required

yes  
Power Processor II/—/yes  
300-450/16, 13 × 100; 16, 13 × 75, Sarstedt  
no  
yes  
yes  
weekly

Automated input/accessioning available

- Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput\*\*
- Containers device accommodates/Dedicated lanes for stat samples
- Maximum No. of samples that can be loaded/Maintenance required

yes  
Power Processor II/—/yes/900  
16, 13 × 100; 16, 13 × 75, Sarstedt/yes  
200/monthly

Automated decapping available

- Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput\*\*
- Containers device accommodates/Maintenance required
- Removes multiple size tube caps per run/Removes screw type sample caps

yes  
Power Processor II/—/yes/600  
16, 13 × 100; 16, 13 × 75, Sarstedt/monthly  
yes/no

Automated sorting available

- Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput\*\*
- Containers device accommodates/Software can sort by

yes  
Power Processor II/—/yes/500  
16, 13 × 100; 16, 13 × 75, Sarstedt/method, output

Specimen integrity monitor available

- Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput\*\*
- Containers device accommodates/Maintenance required

yes  
Power Processor II/—/yes/90  
16, 13 × 100; 16, 13 × 75, Sarstedt/monthly

Automated aliquotting available

- Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput\*\*
- Containers device accommodates
- Inspects samples for bar code/Detects and reports clots in specimen
- Detects and reports quantity not sufficient specimens/Maintenance required

yes  
Power Processor II/—/yes/140 primary samples  
16, 13 × 100; 16, 13 × 75, Sarstedt  
yes/yes  
yes/daily, weekly

Instrument (analyzer) interfaces

- Rules-based instrument interface control subsystem
- Process control of instrument via control subsystem
- Physical/hardware (instrument/specimen) interface
- Hematology/Chemistry/Coagulation

yes  
yes  
robotic arm interface/point-of-reference sampling, robotic arm interface/  
point-of-reference sampling, robotic arm interface  
point-of-reference sampling, robotic arm interface/point-of-reference sampling

- Immunoassay/Urinalysis

Instruments to which your system or product is interfaced

Abbott Architect, Axsym; Siemens Advia, Atlas; Beckman Coulter LX 20, DxC, Dxl; Ortho 950, 250, Eci; Roche Modular; Stago Star

Other robotic products/components to which system or product is linked

—

Automated recapper or sealer available

- Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput\*
- Recaps-seals multiple size tubes simult./Containers device accommodates
- Maintenance required

yes  
Power Processor III/—/yes/500  
no/13 × 100; 13 × 75, Sarstedt  
weekly

Automated storage and retrieval available

- Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput\*
- Containers device accommodates/Connects to the track
- Room temperature/Minimum and maximum number of tubes stored per module
- Multiple size tubes can be stored in the same module/Maintenance required
- Refrigerated storage and retrieval capability
- Longitudinal upgrade pathway or plan to protect users' investments
- Average time to install/Who provides service, support/Hours support is available
- On-site biomedical engineer required/User group meets regularly

yes  
Power Processor III/—/yes/500  
13 × 100; 13 × 75, Sarstedt/yes  
yes/1 and 6,000  
no/weekly  
yes  
Power Processor is expandable for upgrades as lab needs grow  
7-21 days/Beckman Coulter/24-7  
no/yes

List price

- Individual list prices for components
- Process control SW/Transportation systems/Auto. centrifugation
- Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval
- Specimen integrity monitor/Automated aliquot
- Instrument (analyzer) interfaces/Automated recap

varies

Distinguishing features (supplied by company)

\* For basic building block unit

\*\* Average throughput in specimen containers per hour per device

Note: a dash in lieu of an answer means company did not answer question or question is not applicable

refrigerated storage with recapping and auto rerun; totally open system;  
intelligent aliquotting; proven consistent turnaround time results

## Laboratory automation systems and workcells

<p><b>Part 5 of 14</b></p> <p>See <a href="http://captodayonline.com/productguides">captodayonline.com/productguides</a> for an interactive version of guide</p>	<p><b>LABOTIX Automation</b>                  Peter J. Manes peter.manes@labotix.com                  2323 S. 171st Street, Omaha, NE 68130                  402-398-2274 www.labotix.com</p>	<p><b>m-u-t America</b>                  Karsten Wittmann kwittmann@mut-group.com                  3931 Deep Rock Road, Henrico, VA 23233                  804-620-4029 www.mut-group.com</p>
<p><b>Name of system/First year installed/No. of 2011 contracts signed</b>  <b>No. of live sites installed in N. America/Europe/Asia-Australia</b></p>	<p>RRUSH/1994/1                  11/4/0</p>	<p>HCTS2000 MK3 racking device/2008/—                  —</p>
<p><b>Automation products that are available</b></p> <ul style="list-style-type: none"> <li>• Pre-analytical processor/Total laboratory automation</li> <li>• Automated functions: Accessioning/Track load/Centrifugation/Decapping</li> <li>• Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing</li> <li>• Automated functions: Storage-retrieval/Intelligent sample routing</li> <li>• SW: Dedicated Process Control/Middleware control using LIS/Architecture</li> <li>• Company has dedicated automation support team/Remote system monitoring</li> </ul>	<p>yes/yes                  yes/yes/yes/yes                  yes/yes/yes/yes (recapping)                  yes/yes                  yes/yes/open                  yes/yes</p>	<p>yes/no                  yes/no/no/no                  yes/no/no/no                  no/yes                  yes/yes/closed                  yes/yes</p>
<p><b>Software features/functionality</b></p> <ul style="list-style-type: none"> <li>• Patient demographics and insurance data/Rules-based architecture</li> <li>• Supports data retrieval/Internet connectivity</li> <li>• Online real-time help system/QC/Stats and management reports</li> <li>• Evaluates validity and releasability of results from automated analyzers</li> <li>• Specimen tracking/Priority processing/Random-access spec. movement</li> <li>• Supports accession number redundancy (duplicate specimen ID)</li> <li>• Supports specimen carrier and level identification</li> <li>• Unique bar-code number per container required</li> <li>• Specimen routing/Multistop routing (one tube to multiple workstations)</li> <li>• Specimen scheduling/Instrument scheduling</li> <li>• Routes test to workstation/Automatic reflex, repeat, dilutions</li> <li>• Supports multiple HW configuration/Supports other proprietary transport. HW</li> <li>• Sample storage and retrieval SW/Supports approved CLSI standards</li> </ul>	<p>—/automation SW feature                  automation SW feature/—                  automation SW feature/automation SW feature/automation SW feature                  LIS feature                  automation SW feature/automation SW feature/automation SW feature                  LIS feature                  —                  automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/LIS feature                  automation SW feature/LIS feature                  automation SW feature/—                  automation SW feature/automation SW feature</p>	<p>LIS feature/automation SW feature                  —                  automation SW feature/—/—                  —                  automation SW feature/automation SW feature/—                  automation SW feature                  —                  —                  automation SW feature/automation SW feature                  —                  —                  automation SW feature/—                  —/automation SW feature</p>
<p><b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b></p>	<p>Cerner, PGP, Triple G, Sunquest, Rubicon/HL7 or ASTM</p>	<p>Mysis, Soft, DI, VA, DHCP/ASTM</p>
<p><b>Transportation systems available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Containers device accommodates/Average throughput in cm per second</li> <li>• Supports automatic rerouting for reflex-repeat-dilutions</li> <li>• Modular HW/Installed options/Device can operate in track and manual mode</li> <li>• Required utilities/Required maintenance</li> <li>• Carrier type/Scalable system</li> </ul>	<p>yes                  Flexlink/custom by site/yes                  16, 13 × 100; 16, 13 × 75/variable                  yes                  yes/floor mounted, overhead mounted/—                  electricity/quarterly                  single specimen container per carrier/yes</p>	<p>no                  —                  —                  —                  —                  —</p>
<p><b>Automated centrifugation available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Maximum throughput/Containers device accommodates</li> <li>• Can identify tube types for custom programmed rate and spin times per run</li> <li>• More than one centrifuge can be connected to track system</li> <li>• For multi-unit centrifuge, each centrifuge operates independently for rate and time</li> <li>• Maintenance required</li> </ul> <p><b>Automated input/accessioning available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Dedicated lanes for stat samples</li> <li>• Maximum No. of samples that can be loaded/Maintenance required</li> </ul> <p><b>Automated decapping available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> <li>• Removes multiple size tube caps per run/Removes screw type sample caps</li> </ul> <p><b>Automated sorting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Software can sort by</li> </ul> <p><b>Specimen integrity monitor available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> </ul> <p><b>Automated aliquotting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates</li> <li>• Inspects samples for bar code/Detects and reports clots in specimen</li> <li>• Detects and reports quantity not sufficient specimens/Maintenance required</li> </ul>	<p>yes                  Hettich/74 × 34 × 36 inches/yes                  350/16, 13 × 100; 16, 13 × 75                  no                  yes                  yes                  quarterly                  yes                  Labotix/74 × 34 × 36 inches/yes/1,200                  16, 13 × 100; 16, 13 × 75/yes                  1,200/quarterly                  yes                  Labotix/20 × 9 × 12 inches/yes/400                  16, 13 × 100; 16, 13 × 75/quarterly                  yes/no                  yes                  Labotix/74 × 34 × 36 inches/yes/400                  16, 13 × 100; 16, 13 × 75/specimen, method, output                  yes                  —                  —                  yes                  Labotix/60 × 57 × 25 inches/yes/300                  16, 13 × 100; 16, 13 × 75                  yes/yes                  yes/quarterly</p>	<p>no                  —                  —                  —                  —                  —                  yes                  HCTS2000 MK3/61 × 98 × 53 inches/yes/800–2,000                  16, 13 × 100; 16, 13 × 75, 8–19 mm diameter × 75–120 mm height/no                  550/daily, monthly                  no                  —                  —                  —                  yes                  HCTS2000 MK3/61 × 98 × 53 inches/yes/800–2,000                  16, 13 × 100; 16, 13 × 75, 8–19 mm diameter × 75–120 mm height/                  specimen, method, output                  no                  —                  —                  no                  —                  —                  —</p>
<p><b>Instrument (analyzer) interfaces</b></p> <ul style="list-style-type: none"> <li>• Rules-based instrument interface control subsystem</li> <li>• Process control of instrument via control subsystem</li> </ul> <p><b>Physical/hardware (instrument/specimen) interface</b></p> <ul style="list-style-type: none"> <li>• Hematology/Chemistry/Coagulation</li> <li>• Immunoassay/Urinalysis</li> </ul>	<p>yes                  —                  —                  point-of-ref., robotic rack/point-of-ref., robotic rack/point-of-ref., robotic rack                  point-of-reference, robotic rack/point-of-reference, robotic rack</p>	<p>no                  no                  no/no/no                  no/no</p>
<p><b>Instruments to which your system or product is interfaced</b></p> <p><b>Other robotic products/components to which system or product is linked</b></p>	<p>Beckman Coulter Dxl 800, Stago Star Evolution, Olympus 2700 and 5400, Siemens Advia Centaur, Sysmex HST with SMS, Ortho-Clinical Vitros, and more</p> <p>—</p>	<p>—</p> <p>—</p>
<p><b>Automated recapper or sealer available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Recaps-seals multiple size tubes simul./Containers device accommodates</li> <li>• Maintenance required</li> </ul>	<p>recapper                  Labotix/60 × 13 × 23 inches/yes/750                  yes/16, 13 × 100; 16, 13 × 75                  quarterly</p>	<p>no                  —                  —                  —</p>
<p><b>Automated storage and retrieval available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Connects to the track</li> <li>• Room temperature/Minimum and maximum number of tubes stored per module</li> <li>• Multiple size tubes can be stored in the same module/Maintenance required</li> <li>• Refrigerated storage and retrieval capability</li> </ul> <p><b>Longitudinal upgrade pathway or plan to protect users' investments</b></p> <p><b>Average time to install/Who provides service, support/Hours support is available</b></p> <p><b>On-site biomedical engineer required/User group meets regularly</b></p>	<p>yes                  Labotix/90 × 47 × 56 inches/yes/750                  16, 13 × 100; 16, 13 × 75/yes                  no/5,700                  yes/quarterly                  yes                  open system allows changing analyzers and vendors at any time; expands and changes physical shape of track at any time                  30 days/Labotix/24–7–365 days per year                  —/no</p>	<p>no                  —                  —                  —                  —                  independent of analyzer company; module can be upgraded with options                  &lt;2 days/m-u-t America/24-7                  no/no</p>
<p><b>List price</b></p> <p><b>Individual list prices for components</b></p> <ul style="list-style-type: none"> <li>• Process control SW/Transportation systems/Auto. centrifugation</li> <li>• Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval</li> <li>• Specimen integrity monitor/Automated aliquot</li> <li>• Instrument (analyzer) interfaces/Automated recap</li> </ul>	<p>varies                  —                  —                  —                  —</p>	<p>\$161,600                  included/—/—                  —/—/included/—                  —                  —</p>
<p><b>Distinguishing features (supplied by company)</b></p> <p><i>* For basic building block unit</i></p> <p><i>** Average throughput in specimen containers per hour per device</i></p> <p><i>Note: a dash in lieu of an answer means company did not answer question or question is not applicable</i></p>	<p>open system sorts and delivers all specimens to all vendors' analyzers anywhere in lab; users can change vendors without changing automation; scalable system allows customers to expand and reconfigure automation at any time</p>	<p>bulk loading of tubes; tubes are placed into analyzer racks; sorting to output bins and analyzer racks</p>

Tabulation does not represent an endorsement by the College of American Pathologists.

## Laboratory automation systems and workcells

<b>Part 6 of 14</b>	<b>m-u-t America</b> Karsten Wittmann kwittmann@mut-group.com 3931 Deep Rock Road, Henrico, VA 23233 804-620-4029 www.mut-group.com	<b>Ortho-Clinical Diagnostics</b> Mark Steelman msteelma@its.jnj.com 1001 US Route 202, Raritan, NJ 08869 585-453-3420 www.orthoclinical.com
See captodayonline.com/productguides for an interactive version of guide		
<b>Name of system/First year installed/No. of 2011 contracts signed No. of live sites installed in N. America/Europe/Asia-Australia</b>	HCTS2000 MK2 automated sorter/2007/— —	enGen Laboratory Automation System/2001/18 25/71/2
<b>Automation products that are available</b> • Pre-analytical processor/Total laboratory automation • Automated functions: Accessioning/Track load/Centrifugation/Decapping • Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing • Automated functions: Storage-retrieval/Intelligent sample routing • SW: Dedicated Process Control/Middleware control using LIS/Architecture • Company has dedicated automation support team/Remote system monitoring	yes/no yes/no/no/no yes/no/no/no no/yes yes/yes/closed yes/yes	yes/yes yes/yes/yes/yes yes/yes/no/— in development/yes yes/yes/open yes/yes
<b>Software features/functionality</b> • Patient demographics and insurance data/Rules-based architecture • Supports data retrieval/Internet connectivity • Online real-time help system/QC/Stats and management reports • Evaluates validity and releasability of results from automated analyzers • Specimen tracking/Priority processing/Random-access spec. movement • Supports accession number redundancy (duplicate specimen ID) • Supports specimen carrier and level identification • Unique bar-code number per container required • Specimen routing/Multistop routing (one tube to multiple workstations) • Specimen scheduling/Instrument scheduling • Routes test to workstation/Automatic reflex, repeat, dilutions • Supports multiple HW configuration/Supports other proprietary transport. HW • Sample storage and retrieval SW/Supports approved CLSI standards	LIS feature/automation SW feature — automation SW feature/—/— — automation SW feature/automation SW feature/— automation SW feature — — automation SW feature/automation SW feature — — automation SW feature/— —/automation SW feature	automation SW feature/automation SW feature automation SW feature/automation SW feature —/automation SW feature/automation SW feature automation SW feature automation SW feature/automation SW feature/automation SW feature — automation SW feature automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/—
<b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b>	Mysis, Soft, DI, VA, DHCP/ASTM	Cerner, Misys, SCC, several others/HL7, ASTM
<b>Transportation systems available</b> • Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Containers device accommodates/Average throughput in cm per second • Supports automatic rerouting for reflex-repeat-dilutions • Modular HW/Installed options/Device can operate in track and manual mode • Required utilities/Required maintenance • Carrier type/Scalable system	no — — — — —	yes Covered Conveyor/600–2,400 mm sections/yes 16, 13 × 100; 16, 13 × 75/10 yes yes/floor mounted/yes compressed air, electricity/annually single specimen container per carrier/yes
<b>Automated centrifugation available</b> • Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Maximum throughput/Containers device accommodates • Can identify tube types for custom programmed rate and spin times per run • More than one centrifuge can be connected to track system • For multi-unit centrifuge, each centrifuge operates independently for rate and time • Maintenance required	no — — — — —	yes centrifuge module/1,900 × 1,200 × 1,375 mm/yes 400; 96-tube capacity/13 × 100; 13 × 75 yes yes yes quarterly
<b>Automated input/accessioning available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Dedicated lanes for stat samples • Maximum No. of samples that can be loaded/Maintenance required	yes HTS2000 MK2/48 × 56 × 31 inches/yes/2,000 16, 13 × 100; 16, 13 × 75, 8–19 mm diameter × 75–120 mm height/no 550/daily, monthly	yes rack entry-exit module/1,900 × 1,200 × 965 mm/yes/500 16, 13 × 100; 16, 13 × 75/yes 600/annually
<b>Automated decapping available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required • Removes multiple size tube caps per run/Removes screw type sample caps	no — — —	yes decapper module/1,600 × 600 × 965 mm/yes/600 16, 13 × 100; 16, 13 × 75/annually yes/yes
<b>Automated sorting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Software can sort by	yes HCTS2000 MK2/48 × 56 × 31 inches/yes/2,000 16, 13 × 100; 16, 1 × 75, others/specimen type, bar code, cap color, method, others	yes rack exit-entry module/1,900 × 1,200 × 965 mm/yes/500 16, 13 × 100; 16, 13 × 75/specimen, method, output yes via Vitros 5,1 FS 3600, 5600/—/—/— 16, 13 × 100; 16, 13 × 75/weekly, monthly, annually
<b>Specimen integrity monitor available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required	no — —	yes aliquoter and labeler module/1,900 × 1,500 × 965 mm/yes/200 16, 13 × 100; 16, 13 × 75 yes/yes yes/quarterly
<b>Automated aliquotting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates • Inspects samples for bar code/Detects and reports clots in specimen • Detects and reports quantity not sufficient specimens/Maintenance required	no — — — —	yes aliquoter and labeler module/1,900 × 1,500 × 965 mm/yes/200 16, 13 × 100; 16, 13 × 75 yes/yes yes/quarterly
<b>Instrument (analyzer) interfaces</b> • Rules-based instrument interface control subsystem • Process control of instrument via control subsystem Physical/hardware (instrument/specimen) interface • Hematology/Chemistry/Coagulation • Immunoassay/Urinalysis	no no no/no/no no/no	yes — robotic arm interface/point-of-reference sampling/robotic arm interface point-of-reference sampling/—
<b>Instruments to which your system or product is interfaced</b>	—	Vitros 5600, 4600, 3600, 5,1 FS systems; interfaces with coagulation and hematology systems
<b>Other robotic products/components to which system or product is linked</b>	—	—
<b>Automated recapper or sealer available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput* • Recaps-seals multiple size tubes simult./Containers device accommodates • Maintenance required	no — — —	recapper recapper module/1,600 × 600 × 965 mm/yes/500 yes/16, 13 × 100; 16, 13 × 75 annually
<b>Automated storage and retrieval available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput* • Containers device accommodates/Connects to the track • Room temperature/Minimum and maximum number of tubes stored per module • Multiple size tubes can be stored in the same module/Maintenance required • Refrigerated storage and retrieval capability Longitudinal upgrade pathway or plan to protect users' investments	no — — — — — independent of analyzer company; module can be upgraded with options	yes, in development — — — — — customized automation offering, enGen can be reconfigured or upgraded as needs change; SW configuration updates available periodically depends on configurable customizations/depends on service contract with Ortho no/no
<b>Average time to install/Who provides service, support/Hours support is available On-site biomedical engineer required/User group meets regularly</b>	<2 days/m-u-t America/24–7 no/no	
<b>List price</b> Individual list prices for components • Process control SW/Transportation systems/Auto. centrifugation • Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval • Specimen integrity monitor/Automated aliquot • Instrument (analyzer) interfaces/Automated recap	\$116,000 included—/— —/—/included/— — —	varies — — — —
<b>Distinguishing features (supplied by company)</b> * For basic building block unit ** Average throughput in specimen containers per hour per device Note: a dash in lieu of an answer means company did not answer question or question is not applicable	no robotic arms used, high-throughput yields and reliability with ease of operation and installation; users can pour sample tubes into hopper, eliminating shuffling of tubes in and out of racks in lab reception areas; simplicity and flexibility of sorting rules and methods	customizable: systems designed to fit in existing floor space while providing Lean workflow; configurable: systems designed to interface with several lab analyzers; systems grow with the lab



## Laboratory automation systems and workcells

<b>Part 7 of 14</b>	<b>Roche Diagnostics</b> Ed Duning ed.duning@roche.com 9115 Hauge Drive, Indianapolis, IN 46250 317-521-4710 www.roche-diagnostics.us	<b>Roche Diagnostics</b> Ed Duning ed.duning@roche.com 9115 Hauge Drive, Indianapolis, IN 46250 317-521-4710 www.roche-diagnostics.us
See captodayonline.com/productguides for an interactive version of guide		
<b>Name of system/First year installed/No. of 2011 contracts signed No. of live sites installed in N. America/Europe/Asia-Australia</b>	<b>Aliquoting System cobas p612/2002/15 42/165/59</b>	<b>Workstation cobas p612 and cobas p512 connected to EC1/2003/5 5/25/3</b>
<b>Automation products that are available</b> • Pre-analytical processor/Total laboratory automation • Automated functions: Accessioning/Track load/Centrifugation/Decapping • Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing • Automated functions: Storage-retrieval/Intelligent sample routing • SW: Dedicated Process Control/Middleware control using LIS/Architecture • Company has dedicated automation support team/Remote system monitoring	yes/yes yes/yes (as option)/yes/yes yes/yes/yes/yes yes/yes yes/yes/closed yes/yes	yes/yes yes/yes/yes/yes yes/yes/yes/yes yes/yes yes/yes/closed yes/yes
<b>Software features/functionality</b> • Patient demographics and insurance data/Rules-based architecture • Supports data retrieval/Internet connectivity • Online real-time help system/QC/Stats and management reports • Evaluates validity and releasability of results from automated analyzers • Specimen tracking/Priority processing/Random-access spec. movement • Supports accession number redundancy (duplicate specimen ID) • Supports specimen carrier and level identification • Unique bar-code number per container required • Specimen routing/Multistop routing (one tube to multiple workstations) • Specimen scheduling/Instrument scheduling • Routes test to workstation/Automatic reflex, repeat, dilutions • Supports multiple HW configuration/Supports other proprietary transport. HW • Sample storage and retrieval SW/Supports approved CLSI standards	automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature/automation SW feature — automation SW feature/automation SW feature/automation SW feature automation SW feature automation SW feature — automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/— automation SW feature/automation SW feature automation SW feature/automation SW feature	automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature/automation SW feature — automation SW feature/automation SW feature/automation SW feature automation SW feature automation SW feature — automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/— automation SW feature/automation SW feature automation SW feature/automation SW feature
<b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b>	Cerner, MCS, Medat, Systek, MIPS, Providens, Bayer, Molis, Omega, Misys, Vertex, Zanacore, DI, Cirrus, SCC Soft, Nyantech, MCS Promed, Swisslab, Melos, IDAA, Syscomp, OSM, Star LIMS, others/ASTM and system-specific dynamic interface	Cerner, MCS, Medat, Systek, MIPS, Providens, Bayer, Molis, Omega, Misys, Vertex, Zanacore, DI, Cirrus, SCC Soft, Nyantech, MCS Promed, Swisslab, Melos, IDAA, Syscomp, OSM, Star LIMS, others/ASTM and system-specific dynamic interface
<b>Transportation systems available</b> • Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Containers device accommodates/Average throughput in cm per second • Supports automatic rerouting for reflex-repeat-dilutions • Modular HW/Installed options/Device can operate in track and manual mode • Required utilities/Required maintenance • Carrier type/Scalable system	yes transport built into the instrument/—/yes 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 mm up to 15.5 × 108 mm/— no yes/—/yes electricity/weekly, quarterly single specimen container per carrier/yes	yes transport built into the instrument/—/yes 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 mm up to 15.5 × 108 mm/— no yes/floor mounted/yes compressed air, electricity/weekly, quarterly single and multiple (5) specimen container per carrier/yes
<b>Automated centrifugation available</b> • Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5  • Maximum throughput/Containers device accommodates • Can identify tube types for custom programmed rate and spin times per run • More than one centrifuge can be connected to track system • For multi-unit centrifuge, each centrifuge operates independently for rate and time • Maintenance required	yes single (EC1)/61.4 × 78.3 × 83.6 inches; EC2: 85.8 × 79.3 × 78.7 inches/yes  EC1: 380 tubes per hour/16, 13 × 100; 16, 13 × 75, others yes yes yes weekly, quarterly	yes single (EC1)/61.4 × 78.3 × 83.6 inches/yes  EC1: 380 tubes per hour/16, 13 × 100; 16, 13 × 75, others yes yes yes daily, quarterly
<b>Automated input/accessioning available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Dedicated lanes for stat samples • Maximum No. of samples that can be loaded/Maintenance required	yes input unit as part of system/78.74 × 33.47 × 69.29 inches/yes/up to 1,200 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 mm up to 15.5 × 108 mm/yes 600/daily, quarterly	yes input unit as part of instrument/78.74 × 33.47 × 69.29 inches/yes/up to 1,200 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 mm up to 15.5 × 108 mm/yes EC1: 150 tubes/daily, quarterly
<b>Automated decapping available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required • Removes multiple size tube caps per run/Removes screw type sample caps	yes decapping module as part of system/14.96 × 12.60 × 5.90 inches/yes/up to 1,200 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm/daily, quarterly yes/yes	yes decapping module as part of instrument/14.96 × 12.60 × 5.90 in./yes/up to 1,200 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm/daily, quarterly yes/yes
<b>Automated sorting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Software can sort by	yes output sorter as part of system/71.65 × 55.90 × 55.11 inches/yes/up to 1,200 16, 13×100; 16, 13×75; 11.5×65.5 to 15.5×108 mm/specimen, method, output	yes part of Aliquoting System cobas p612 or Sorting System cobas p512/—/yes/up to 1,200 16, 13×100; 16, 13×75; 11.5×65.5 to 15.5×108 mm/specimen, method, output
<b>Specimen integrity monitor available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required	yes Quality Check Unit QS I/11.4 × 19.7 × 14.0 inches/yes/850 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm/daily, quarterly	yes Quality Check Unit QS I/11.4 × 19.7 × 14.0 inches/yes/850 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm/daily, quarterly
<b>Automated aliquotting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates • Inspects samples for bar code/Detects and reports clots in specimen • Detects and reports quantity not sufficient specimens/Maintenance required	yes aliquoting unit as part of system/125 × 73.2 × 78.7 inches/yes/655 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm yes/yes yes/daily, quarterly	yes aliquoting unit as part of Aliquoting System cobas p612/125 × 73.2 × 78.7 in./yes/655 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm yes/yes yes/daily, quarterly
<b>Instrument (analyzer) interfaces</b> • Rules-based instrument interface control subsystem • Process control of instrument via control subsystem Physical/hardware (instrument/specimen) interface • Hematology/Chemistry/Coagulation • Immunoassay/Urinalysis	yes no  no/no/no no/no	yes no  no/no/no no/no
<b>Instruments to which your system or product is interfaced</b>	—	—
<b>Other robotic products/components to which system or product is linked</b>	—	—
<b>Automated recapper or sealer available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Recaps-seals multiple size tubes simult./Containers device accommodates • Maintenance required	sealer recapping module as part of system/13.39 × 12.20 × 8.66 inches/yes/up to 1,200 yes/16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm daily, quarterly	sealer recapping module/13.39 × 12.20 × 8.66 inches/yes/up to 1,200 yes/16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm daily, quarterly
<b>Automated storage and retrieval available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Connects to the track • Room temperature/Minimum and maximum number of tubes stored per module • Multiple size tubes can be stored in the same module/Maintenance required • Refrigerated storage and retrieval capability Longitudinal upgrade pathway or plan to protect users' investments	yes as part of system (output sorter), up to 41 workplaces/—/yes/up to 1,200 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm/yes no/1,200 yes/daily, quarterly no independent of any analyzer company, Roche/PVT modules can be upgraded	yes implemented into system, up to 41 workplaces/—/yes/up to 1,200 16, 13 × 100; 16, 13 × 75; 11.5 × 65.5 to 15.5 × 108 mm/no no/1,200 yes/— no independent of any analyzer company, Roche/PVT modules can be upgraded
<b>Average time to install/Who provides service, support/Hours support is available On-site biomedical engineer required/User group meets regularly</b>	~1-2 weeks/Roche Diagnostics/daily 8 AM-5 PM (EST) and 24-7 upon request no/no	~1-2 weeks/Roche Diagnostics/daily 8 AM-5 PM (EST) and 24-7 upon request no/no
<b>List price</b> Individual list prices for components • Process control SW/Transportation systems/Auto. centrifugation • Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval • Specimen integrity monitor/Automated aliquot • Instrument (analyzer) interfaces/Automated recap	— —/included/—/included included/included/included/— —/included —	— —/included/included/included included/included/included/— —/included as part of Aliquoting System cobas p612 —
<b>Distinguishing features (supplied by company)</b> * For basic building block unit ** Average throughput in specimen containers per hour per device	basic platform can be configured for each customer routine workflow using many vendor sample carriers for input and output sorting and archiving; recursive workflow allows samples to be processed multiple times	basic platform can be configured for each customer routine workflow using many vendor sample carriers for input and output sorting and archiving; recursive workflow allows samples to be processed multiple times

## Laboratory automation systems and workcells

<p><b>Part 8 of 14</b></p> <p>See <a href="http://captodayonline.com/productguides">captodayonline.com/productguides</a> for an interactive version of guide</p>	<p><b>Roche Diagnostics</b>                  Ed Duning ed.duning@roche.com                  9115 Hauge Drive, Indianapolis, IN 46250                  317-521-4710 www.roche-diagnostics.us</p>	<p><b>Roche Diagnostics</b>                  Ed Duning ed.duning@roche.com                  9115 Hauge Drive, Indianapolis, IN 46250                  317-521-4710 www.roche-diagnostics.us</p>
<p><b>Name of system/First year installed/No. of 2011 contracts signed</b>  <b>No. of live sites installed in N. America/Europe/Asia-Australia</b></p>	<p><b>Sorting System cobas p512/2001/23</b>                  28/112/20</p>	<p><b>Modular Pre-Analytics EVO/2000/72</b>                  172/353/265</p>
<p><b>Automation products that are available</b></p> <ul style="list-style-type: none"> <li>• Pre-analytical processor/Total laboratory automation</li> <li>• Automated functions: Accessioning/Track load/Centrifugation/Decapping</li> <li>• Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing</li> <li>• Automated functions: Storage-retrieval/Intelligent sample routing</li> <li>• SW: Dedicated Process Control/Middleware control using LIS/Architecture</li> <li>• Company has dedicated automation support team/Remote system monitoring</li> </ul>	<p>yes/yes                  yes/yes/yes/yes                  yes/no/no/yes                  yes/yes                  yes/yes/closed                  yes/yes</p>	<p>yes/yes                  yes/yes/yes/yes                  yes/yes/yes/yes                  yes/yes                  yes/yes/open and closed                  yes/yes</p>
<p><b>Software features/functionality</b></p> <ul style="list-style-type: none"> <li>• Patient demographics and insurance data/Rules-based architecture</li> <li>• Supports data retrieval/Internet connectivity</li> <li>• Online real-time help system/QC/Stats and management reports</li> <li>• Evaluates validity and releasability of results from automated analyzers</li> <li>• Specimen tracking/Priority processing/Random-access spec. movement</li> <li>• Supports accession number redundancy (duplicate specimen ID)</li> <li>• Supports specimen carrier and level identification</li> <li>• Unique bar-code number per container required</li> <li>• Specimen routing/Multistop routing (one tube to multiple workstations)</li> <li>• Specimen scheduling/Instrument scheduling</li> <li>• Routes test to workstation/Automatic reflex, repeat, dilutions</li> <li>• Supports multiple HW configuration/Supports other proprietary transport. HW</li> <li>• Sample storage and retrieval SW/Supports approved CLSI standards</li> </ul>	<p>automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature/automation SW feature                  —                  automation SW feature/automation SW feature/automation SW feature                  automation SW feature                  automation SW feature                  —                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/—                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature</p>	<p>automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature/automation SW feature                  automation SW feature                  automation SW feature/automation SW feature/automation SW feature                  automation SW feature                  automation SW feature                  automation SW feature                  automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/—                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature</p>
<p><b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b></p>	<p>Cerner, MCS, Medat, Syspek, MIPS, Providens, Bayer, Molis, Omega, Misys, Vertex, Zanacore, DI, Cirrus, SCC Soft, Nyantech, MCS Promed, Swisslab, Melos, IDAA, Syscomp, OSM, Star LIMS, others/ASTM and system-specific dynamic interface</p>	<p>Cerner, MCS, Medat, Syspek, MIPS, Providens, Bayer, Molis, Omega, Misys, Vertex, Zanacore, DI, Cirrus, SCC Soft, Nyantech, MCS Promed, Swisslab, Melos, IDAA, Syscomp, OSM, Star LIMS, others/LIS to LAS, HL7, ASTM</p>
<p><b>Transportation systems available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions* (H x W x D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Containers device accommodates/Average throughput in cm per second</li> <li>• Supports automatic rerouting for reflex-repeat-dilutions</li> <li>• Modular HW/Installed options/Device can operate in track and manual mode</li> <li>• Required utilities/Required maintenance</li> <li>• Carrier type/Scalable system</li> </ul>	<p>yes                  transport built into the instrument/—/yes                  16, 13 x 100; 16, 13 x 75; 11.5 x 65.5 to 15.5 x 108 mm/—                  no                  yes/—/yes                  electricity/daily, quarterly                  single specimen container per carrier/yes</p>	<p>yes                  MPA (A, B, C)/A: 4.6 x 15 x 3.5 ft.; B: 4.6 x 18 x 3.5 ft.; C: 4.6 x 9 x 3.5 feet/yes                  16, 13 x 100; 16, 13 x 75; 13 x 92, Greiner FBT, others/400 tubes per hour                  no                  yes/floor mounted/yes                  electricity/daily, quarterly                  multiple specimen (5) container per carrier/yes</p>
<p><b>Automated centrifugation available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions (H x W x D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Maximum throughput/Containers device accommodates</li> <li>• Can identify tube types for custom programmed rate and spin times per run</li> <li>• More than one centrifuge can be connected to track system</li> <li>• For multi-unit centrifuge, each centrifuge operates independently for rate and time</li> <li>• Maintenance required</li> </ul> <p><b>Automated input/accessioning available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Dedicated lanes for stat samples</li> <li>• Maximum No. of samples that can be loaded/Maintenance required</li> </ul> <p><b>Automated decapping available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> <li>• Removes multiple size tube caps per run/Removes screw type sample caps</li> </ul> <p><b>Automated sorting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Software can sort by</li> </ul> <p><b>Specimen integrity monitor available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> </ul> <p><b>Automated aliquotting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates</li> <li>• Inspects samples for bar code/Detects and reports clots in specimen</li> <li>• Detects and reports quantity not sufficient specimens/Maintenance required</li> </ul>	<p>yes                  single (EC1)/61.4x78.3x83.6 inches; EC2: 85.8x79.3x78.7 in./yes                  EC1: 380 tubes per hour/16, 13 x 100; 16, 13 x 75, others                  yes                  yes                  yes                  yes                  daily, quarterly                  yes                  input unit as part of system/78.74 x 33.47 x 69.29 inches/yes/up to 1,200                  16, 13 x 100; 16, 13 x 75; 11.5 x 65.5 mm up to 15.5 x 108 mm/yes                  600/daily, quarterly                  yes                  decapping module as part of system/14.96 x 12.60 x 5.90 inches/yes/up to 1,200                  16, 13 x 100; 16, 13 x 75; 11.5 x 65.5 to 15.5 x 108 mm/daily, quarterly                  yes/yes                  yes                  output sorter as part of system/71.65 x 55.90 x 55.11 inches/yes/up to 1,200                  16, 13x100; 16, 13x75; 11.5x65.5 to 15.5x108 mm/specimen, method, output                  yes                  Quality Check Unit QS I/11.4 x 19.7 x 14.0 inches/yes/850                  16, 13 x 100; 16, 13 x 75; 11.5 x 65.5 to 15.5 x 108 mm/daily, quarterly                  no                  —                  —                  —                  —</p>	<p>yes                  standard centrifuge/3 x 2.5 x 3.5 feet/yes                  250/16, 13 x 100; 16, 13 x 75                  yes                  yes                  no                  daily, quarterly                  yes                  standard input buffer/42 x 38 x 41 inches/yes/160 racks                  16, 13 x 100; 16, 13 x 75/yes                  300/daily, quarterly                  yes                  standard decapper/49 x 18 x 41 inches/yes/80 racks                  16, 13x100; 16, 13x75; rubber, hemogard, twist-off/daily, quarterly                  yes/yes                  yes                  standard sorter/36.6 x 11.8 x 41 inches/yes/80 racks                  16, 13 x 100; 16, 13 x 75; 13x92, Greiner FBT, others/specimen, method, output                  yes                  standard aliquoter/53 x 42 x 41 inches/yes/80 racks                  16, 13 x 100; 16, 13 x 75; 13x92, Greiner FBT, Greiver, others/daily, quarterly                  yes                  standard aliquoter/53 x 42 x 41 in.ches/yes/80 racks                  16, 13 x 100; 16, 13 x 75; 13x92, Greiner FBT, others                  yes/yes                  yes/daily, quarterly</p>
<p><b>Instrument (analyzer) interfaces</b></p> <ul style="list-style-type: none"> <li>• Rules-based instrument interface control subsystem</li> <li>• Process control of instrument via control subsystem</li> </ul> <p><b>Physical/hardware (instrument/specimen) interface</b></p> <ul style="list-style-type: none"> <li>• Hematology/Chemistry/Coagulation</li> <li>• Immunoassay/Urinalysis</li> </ul>	<p>yes                  no                  no/no/no                  no/no</p>	<p>yes                  yes                  no/point-of-reference sampling/no                  point-of-reference sampling/point-of-reference sampling</p>
<p><b>Instruments to which your system or product is interfaced</b>  <b>Other robotic products/components to which system or product is linked</b></p>	<p>—                  —</p>	<p>Hitachi, Stago                  Hitachi, Stago</p>
<p><b>Automated recapper or sealer available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Recaps-seals multiple size tubes simult./Containers device accommodates</li> <li>• Maintenance required</li> </ul>	<p>sealer                  recapping module as part of system/13.39 x 12.20 x 8.66 inches/yes/up to 1,200                  yes/16, 13 x 100; 16, 13 x 75; 11.5 x 65.5 to 15.5 x 108 mm                  daily, quarterly</p>	<p>recapper                  standard recapper/50 x 17.5 x 41 inches/yes/80 racks                  yes/16, 13 x 100; 16, 13 x 75, 13 x 92, Greiner FBT, Greiver, others                  daily, quarterly</p>
<p><b>Automated storage and retrieval available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Connects to the track</li> <li>• Room temperature/Minimum and maximum number of tubes stored per module</li> <li>• Multiple size tubes can be stored in the same module/Maintenance required</li> <li>• Refrigerated storage and retrieval capability</li> </ul> <p><b>Longitudinal upgrade pathway or plan to protect users' investments</b>  <b>Average time to install/Who provides service, support/Hours support is available</b>  <b>On-site biomedical engineer required/User group meets regularly</b></p>	<p>yes                  archiving included as part of system (output sorter), up to 41                  workplaces/—/yes/up to 1,200                  16, 13 x 100; 16, 13 x 75; 11.5 x 65.5 mm up to 15.5 x 108 mm/yes                  no/up to 1,200                  yes/daily, quarterly                  no                  independent of any analyzer company, Roche/PVT modules can be upgraded                  ~1 week/Roche Diagnostics/daily 8 AM-5 PM (EST); 24-7 upon request                  no/no</p>	<p>yes                  p501, p701/p501: 5.3 ft. x 14 feet; p701: 5.3 feet x 17.6 feet/yes/80 racks                  16, 13 x 100; 16, 13 x 75, 13 x 92, Greiner FBT, Greiver pour-off tube, others/yes                  no/p501: 13,500; p701: 27,000                  yes/daily, quarterly                  yes                  support for a minimum of 10 years after production                  up to 2 weeks/Roche Diagnostics phone and engineering field support/24-7                  no/no</p>
<p><b>List price</b>  <b>Individual list prices for components</b></p> <ul style="list-style-type: none"> <li>• Process control SW/Transportation systems/Auto. centrifugation</li> <li>• Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval</li> <li>• Specimen integrity monitor/Automated aliquot</li> <li>• Instrument (analyzer) interfaces/Automated recap</li> </ul>	<p>—                  —/included/—/included                  included/included/included/—                  —                  —</p>	<p>—                  included/included/included                  included/included/included/included                  included/included                  included/included</p>
<p><b>Distinguishing features (supplied by company)</b>                  * For basic building block unit                  ** Average throughput in specimen containers per hour per device                  Note: a dash in lieu of an answer means company did not answer question or question is not applicable</p>	<p>basic platform can be configured for each customer routine workflow using many vendor sample carriers for input and output sorting and archiving; recursive workflow allows samples to be processed multiple times; quality module QS I for monitoring specimen integrity and measuring volume</p>	<p>scalable and flexible to fit customer needs and facility space requirements; programmed and personalized to customer workflow requirements; three models can be configured in 100+ standard layouts, connecting up to 12 chemistry/immunochemistry modules</p>



## Laboratory automation systems and workcells

<b>Part 9 of 14</b>	<b>Roche Diagnostics</b> Ed Duning ed.duning@roche.com 9115 Hauge Drive, Indianapolis, IN 46250 317-521-4710 www.roche-diagnostics.us	<b>Roche Diagnostics</b> Ed Duning ed.duning@roche.com 9115 Hauge Drive, Indianapolis, IN 46250 317-521-4710 www.roche-diagnostics.us
See captodayonline.com/productguides for an interactive version of guide		
<b>Name of system/First year installed/No. of 2011 contracts signed No. of live sites installed in N. America/Europe/Asia-Australia</b>	<b>cobas p501 (storage and retrieval)/2009/3 1/11/2</b>	<b>cobas p701 (storage and retrieval)/2009/4 4/8/2</b>
<b>Automation products that are available</b> • Pre-analytical processor/Total laboratory automation • Automated functions: Accessioning/Track load/Centrifugation/Decapping • Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing • Automated functions: Storage-retrieval/Intelligent sample routing • SW: Dedicated Process Control/Middleware control using LIS/Architecture • Company has dedicated automation support team/Remote system monitoring	no/no no/no/no/yes no/no/no/yes yes/no yes/yes/closed yes/yes	no/no no/no/no/yes no/no/no/yes yes/no yes/yes/closed yes/yes
<b>Software features/functionality</b> • Patient demographics and insurance data/Rules-based architecture • Supports data retrieval/Internet connectivity • Online real-time help system/QC/Stats and management reports • Evaluates validity and releasability of results from automated analyzers • Specimen tracking/Priority processing/Random-access spec. movement • Supports accession number redundancy (duplicate specimen ID) • Supports specimen carrier and level identification • Unique bar-code number per container required • Specimen routing/Multistop routing (one tube to multiple workstations) • Specimen scheduling/Instrument scheduling • Routes test to workstation/Automatic reflex, repeat, dilutions • Supports multiple HW configuration/Supports other proprietary transport. HW • Sample storage and retrieval SW/Supports approved CLSI standards	—/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature/automation SW feature — automation SW feature/automation SW feature/— automation SW feature — automation SW feature — automation SW feature/— — — automation SW feature/automation SW feature	—/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature/automation SW feature — automation SW feature/automation SW feature/— automation SW feature — automation SW feature — automation SW feature/— — — automation SW feature/automation SW feature
<b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b>	<b>Cerner, Misys, SCC Soft/HL7</b>	<b>Cerner, Misys, SCC Soft/HL7</b>
<b>Transportation systems available</b> • Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Containers device accommodates/Average throughput in cm per second • Supports automatic rerouting for reflex-repeat-dilutions • Modular HW/Installed options/Device can operate in track and manual mode • Required utilities/Required maintenance • Carrier type/Scalable system	no — — — — —	no — — — — —
<b>Automated centrifugation available</b> • Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Maximum throughput/Containers device accommodates • Can identify tube types for custom programmed rate and spin times per run • More than one centrifuge can be connected to track system • For multi-unit centrifuge, each centrifuge operates independently for rate and time • Maintenance required	no — — — — —	no — — — — —
<b>Automated input/accessioning available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Dedicated lanes for stat samples • Maximum No. of samples that can be loaded/Maintenance required	yes manual and connected to Roche automation/14 × 5.3 × 7.5 feet/yes/400 16, 13 × 100; 16, 13 × 75, 11.5 × 65.5 mm—15.5 × 108 mm/— 300 manual and continuous from MPA/daily, quarterly	yes manual and connected to Roche automation/17.5 × 5.3 × 7.5 feet/yes/400 16, 13 × 100; 16, 13 × 75, 11.5 × 65.5 mm—15.5 × 108 mm/— 300 manual and continuous from MPA/daily, quarterly
<b>Automated decapping available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required • Removes multiple size tube caps per run/Removes screw type sample caps	yes decapper as part of system/—/yes/total system is 400 16, 13 × 100; 16, 13 × 75, 11.5 × 65.5 mm—15.5 × 108 mm/daily, quarterly yes/yes	yes decapper as part of system/—/yes/total system is 400 16, 13 × 100; 16, 13 × 75, 11.5 × 65.5 mm—15.5 × 108 mm/daily, quarterly yes/yes
<b>Automated sorting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Software can sort by	yes sorter as part of system/—/yes/total system is 400 16, 13 × 100; 16, 13 × 75; 13 × 92, Greiner FBT, others/specimen, output	yes sorter as part of system/—/yes/total system is 400 16, 13 × 100; 16, 13 × 75; 13 × 92, Greiner FBT, others/specimen, output
<b>Specimen integrity monitor available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required	no — —	no — —
<b>Automated aliquotting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates • Inspects samples for bar code/Detects and reports clots in specimen • Detects and reports quantity not sufficient specimens/Maintenance required	no — — — —	no — — — —
<b>Instrument (analyzer) interfaces</b> • Rules-based instrument interface control subsystem • Process control of instrument via control subsystem <b>Physical/hardware (instrument/specimen) interface</b> • Hematology/Chemistry/Coagulation • Immunoassay/Urinalysis	yes yes no/no/no no/no	yes yes no/no/no no/no
<b>Instruments to which your system or product is interfaced Other robotic products/components to which system or product is linked</b>	<b>Roche MPA</b> —	<b>Roche MPA</b> —
<b>Automated recapper or sealer available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Recaps-seals multiple size tubes simult./Containers device accommodates • Maintenance required	recapper recapper as part of system/—/yes/total system is 400 yes/16, 13 × 100; 16, 13 × 75, 13 × 92, 11.5 × 65.5 mm—15.5 × 108 mm daily, quarterly	recapper recapper as part of system/—/yes/total system is 400 yes/16, 13 × 100; 16, 13 × 75, 13 × 92, 11.5 × 65.5 mm—15.5 × 108 mm daily, quarterly
<b>Automated storage and retrieval available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Connects to the track • Room temperature/Minimum and maximum number of tubes stored per module • Multiple size tubes can be stored in the same module/Maintenance required • Refrigerated storage and retrieval capability <b>Longitudinal upgrade pathway or plan to protect users' investments</b>	yes cobas p501/14 × 5.3 × 7.5 feet/yes/400 16, 13 × 100; 16, 13 × 75, 13 × 92, 11.5 × 65.5 mm—15.5 × 108 mm/yes no/13,500 yes/daily, quarterly yes support for a minimum of 10 years after production; product upgrades installed as required	yes cobas p501/17.5 × 5.3 × 7.5 feet/yes/400 16, 13 × 100; 16, 13 × 75, 13 × 92, 11.5 × 65.5 mm—15.5 × 108 mm/yes no/27,000 yes/daily, quarterly yes support for a minimum of 10 years after production; product upgrades installed as required
<b>Average time to install/Who provides service, support/Hours support is available On-site biomedical engineer required/User group meets regularly</b>	<b>1 week/Roche Diagnostics/24-7</b> no/no	<b>1 week/Roche Diagnostics/24-7</b> no/no
<b>List price</b> <b>Individual list prices for components</b> • Process control SW/Transportation systems/Auto. centrifugation • Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval • Specimen integrity monitor/Automated aliquot • Instrument (analyzer) interfaces/Automated recap	varies included/—/— included/included/included/based on system options — —/included	varies included/—/— included/included/included/based on system options — —/included
<b>Distinguishing features (supplied by company)</b> * For basic building block unit ** Average throughput in specimen containers per hour per device <i>Note: a dash in lieu of an answer means company did not answer question or question is not applicable</i>	<b>13,500-tube storage capacity with multiple storage durations for 13- and 16-mm tubes; automatically disposes of tubes at the expiration of the selected storage duration; accept tubes for storage from an automatic feed and manual walk up</b>	<b>27,000-tube storage capacity with multiple storage durations for 13- and 16-mm tubes; automatically disposes of tubes at the expiration of the selected storage duration; accept tubes for storage from an automatic feed and manual walk up</b>

## Laboratory automation systems and workcells

<p><b>Part 10 of 14</b></p> <p><i>See captodayonline.com/productguides for an interactive version of guide</i></p>	<p><b>Sarstedt, Inc.</b>                  Peter Rumswinkel, VP/GM sarstedt@bellsouth.net                  P. O. Box 468, Newton, NC 28658                  800-257-5101 www.sarstedt.com</p>	<p><b>Sarstedt, Inc.</b>                  Peter Rumswinkel, VP/GM sarstedt@bellsouth.net                  P. O. Box 468, Newton, NC 28658                  800-257-5101 www.sarstedt.com</p>
<p><b>Name of system/First year installed/No. of 2011 contracts signed</b>  <b>No. of live sites installed in N. America/Europe/Asia-Australia</b></p>	<p>DC/RC 900 Flex/2009/— —</p>	<p>HSS High Speed Sorter 1625/2004/— —</p>
<p><b>Automation products that are available</b></p> <ul style="list-style-type: none"> <li>• Pre-analytical processor/Total laboratory automation</li> <li>• Automated functions: Accessioning/Track load/Centrifugation/Decapping</li> <li>• Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing</li> <li>• Automated functions: Storage-retrieval/Intelligent sample routing</li> <li>• SW: Dedicated Process Control/Middleware control using LIS/Architecture</li> <li>• Company has dedicated automation support team/Remote system monitoring</li> </ul>	<p>yes/no yes/no/no/yes yes/no/no/yes no/yes yes/yes/open yes/yes</p>	<p>yes/no yes/no/no/yes yes/no/no/yes no/yes yes/yes/open yes/yes</p>
<p><b>Software features/functionality</b></p> <ul style="list-style-type: none"> <li>• Patient demographics and insurance data/Rules-based architecture</li> <li>• Supports data retrieval/Internet connectivity</li> <li>• Online real-time help system/QC/Stats and management reports</li> <li>• Evaluates validity and releasability of results from automated analyzers</li> <li>• Specimen tracking/Priority processing/Random-access spec. movement</li> <li>• Supports accession number redundancy (duplicate specimen ID)</li> <li>• Supports specimen carrier and level identification</li> <li>• Unique bar-code number per container required</li> <li>• Specimen routing/Multistop routing (one tube to multiple workstations)</li> <li>• Specimen scheduling/Instrument scheduling</li> <li>• Routes test to workstation/Automatic reflex, repeat, dilutions</li> <li>• Supports multiple HW configuration/Supports other proprietary transport. HW</li> <li>• Sample storage and retrieval SW/Supports approved CLSI standards</li> </ul>	<p>—/automation SW feature automation SW feature/automation SW feature —/—/automation SW feature — automation SW feature/automation SW feature/automation SW feature automation SW feature — — automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/— —/automation SW feature</p>	<p>—/automation SW feature automation SW feature/automation SW feature —/—/automation SW feature — automation SW feature/automation SW feature/automation SW feature automation SW feature automation SW feature — automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/— —/automation SW feature</p>
<p><b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b></p>	<p>—</p>	<p>—</p>
<p><b>Transportation systems available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Containers device accommodates/Average throughput in cm per second</li> <li>• Supports automatic rerouting for reflex-repeat-dilutions</li> <li>• Modular HW/Installed options/Device can operate in track and manual mode</li> <li>• Required utilities/Required maintenance</li> <li>• Carrier type/Scalable system</li> </ul>	<p>no — — — — —</p>	<p>no — — — — —</p>
<p><b>Automated centrifugation available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Maximum throughput/Containers device accommodates</li> <li>• Can identify tube types for custom programmed rate and spin times per run</li> <li>• More than one centrifuge can be connected to track system</li> <li>• For multi-unit centrifuge, each centrifuge operates independently for rate and time</li> <li>• Maintenance required</li> </ul> <p><b>Automated input/accessioning available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Dedicated lanes for stat samples</li> <li>• Maximum No. of samples that can be loaded/Maintenance required</li> </ul> <p><b>Automated decapping available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> <li>• Removes multiple size tube caps per run/Removes screw type sample caps</li> </ul> <p><b>Automated sorting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Software can sort by</li> </ul> <p><b>Specimen integrity monitor available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> </ul> <p><b>Automated aliquotting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates</li> <li>• Inspects samples for bar code/Detects and reports clots in specimen</li> <li>• Detects and reports quantity not sufficient specimens/Maintenance required</li> </ul>	<p>no — — — — — yes —/—/yes/800 16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100/yes 600/daily, annually yes —/—/yes/800 16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100/daily, annually yes/yes yes —/—/yes/800 16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100/specimen, method, output no — — — no — — — no — — — —</p>	<p>no — — — — — yes —/—/yes/1,200 16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100/yes 600/daily, annually yes —/—/yes/1,200 16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100/daily, annually yes/yes yes —/—/yes/1,200 16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100/specimen, method, output yes —/—/yes/700 16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100/daily, annually no — — — —</p>
<p><b>Instrument (analyzer) interfaces</b></p> <ul style="list-style-type: none"> <li>• Rules-based instrument interface control subsystem</li> <li>• Process control of instrument via control subsystem</li> </ul> <p><b>Physical/hardware (instrument/specimen) interface</b></p> <ul style="list-style-type: none"> <li>• Hematology/Chemistry/Coagulation</li> <li>• Immunoassay/Urinalysis</li> </ul>	<p>no no — —</p>	<p>no no — —</p>
<p><b>Instruments to which your system or product is interfaced</b>  <b>Other robotic products/components to which system or product is linked</b></p>	<p>— —</p>	<p>— —</p>
<p><b>Automated recapper or sealer available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Recaps-seals multiple size tubes simult./Containers device accommodates</li> <li>• Maintenance required</li> </ul>	<p>recapper —/—/yes/800 yes/16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100 daily, annually</p>	<p>recapper —/—/yes/1,200 yes/16, 13 × 100; 16, 13 × 75; 13 × 65 to 16 × 100 daily, annually</p>
<p><b>Automated storage and retrieval available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Connects to the track</li> <li>• Room temperature/Minimum and maximum number of tubes stored per module</li> <li>• Multiple size tubes can be stored in the same module/Maintenance required</li> <li>• Refrigerated storage and retrieval capability</li> </ul> <p><b>Longitudinal upgrade pathway or plan to protect users' investments</b></p>	<p>no — — — — systems are upgradable</p>	<p>no — — — — systems are upgradable</p>
<p><b>Average time to install/Who provides service, support/Hours support is available</b>  <b>On-site biomedical engineer required/User group meets regularly</b></p>	<p>3 days/Sarstedt/M-F 8:00 AM–5 PM no/no</p>	<p>2 weeks/Sarstedt/M-F 8:00 AM–5 PM no/no</p>
<p><b>List price</b></p> <ul style="list-style-type: none"> <li>• Individual list prices for components</li> <li>• Process control SW/Transportation systems/Auto. centrifugation</li> <li>• Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval</li> <li>• Specimen integrity monitor/Automated aliquot</li> <li>• Instrument (analyzer) interfaces/Automated recap</li> </ul>	<p>— — — — —</p>	<p>— — — — —</p>
<p><b>Distinguishing features (supplied by company)</b>  <i>* For basic building block unit</i>  <i>** Average throughput in specimen containers per hour per device</i>  <i>Note: a dash in lieu of an answer means company did not answer question or question is not applicable</i></p>	<p>small sorter footprint; maximizes floor space; fills a gap experienced by smaller labs when large automation is too expensive; supports multiple runs for routine and archiving</p>	<p>small footprint requires minimal lab space; modular design enables configuration with only the necessary modules and functions; custom sort target and rules are determined by the user</p>

Tabulation does not represent an endorsement by the College of American Pathologists.

## Laboratory automation systems and workcells

<p><b>Part 11 of 14</b></p> <p>See <a href="http://captodayonline.com/productguides">captodayonline.com/productguides</a> for an interactive version of guide</p>	<p><b>Sarstedt, Inc.</b>                  Peter Rumswinkel, VP/GM sarstedt@bellsouth.net                  P. O. Box 468, Newton, NC 28658                  800-257-5101 www.sarstedt.com</p>	<p><b>Siemens Healthcare Diagnostics</b>                  Sepehr Seyedzadeh                  511 Benedict Avenue, Tarrytown, NY 10591                  914-524-3827 www.usa.siemens.com/diagnostics</p>
<p><b>Name of system/First year installed/No. of 2011 contracts signed</b>  <b>No. of live sites installed in N. America/Europe/Asia-Australia</b></p>	<p><b>Sarstedt PVS/—</b>                  —</p>	<p><b>ADVIA Solutions/1998/—</b>                  &gt;150 U.S./&gt;500 worldwide</p>
<p><b>Automation products that are available</b></p> <ul style="list-style-type: none"> <li>• Pre-analytical processor/Total laboratory automation</li> <li>• Automated functions: Accessioning/Track load/Centrifugation/Decapping</li> <li>• Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing</li> <li>• Automated functions: Storage-retrieval/Intelligent sample routing</li> <li>• SW: Dedicated Process Control/Middleware control using LIS/Architecture</li> <li>• Company has dedicated automation support team/Remote system monitoring</li> </ul>	<p>yes/no                  yes/—/no/yes                  yes/yes/yes/yes                  no/yes                  yes/yes/open                  yes/yes</p>	<p>yes/yes                  yes/yes/yes/yes                  yes/no/no/no/ various partnerships in place                  yes/yes                  yes/yes/—                  yes/yes</p>
<p><b>Software features/functionality</b></p> <ul style="list-style-type: none"> <li>• Patient demographics and insurance data/Rules-based architecture</li> <li>• Supports data retrieval/Internet connectivity</li> <li>• Online real-time help system/QC/Stats and management reports</li> <li>• Evaluates validity and releasability of results from automated analyzers</li> <li>• Specimen tracking/Priority processing/Random-access spec. movement</li> <li>• Supports accession number redundancy (duplicate specimen ID)</li> <li>• Supports specimen carrier and level identification</li> <li>• Unique bar-code number per container required</li> <li>• Specimen routing/Multistop routing (one tube to multiple workstations)</li> <li>• Specimen scheduling/Instrument scheduling</li> <li>• Routes test to workstation/Automatic reflex, repeat, dilutions</li> <li>• Supports multiple HW configuration/Supports other proprietary transport. HW</li> <li>• Sample storage and retrieval SW/Supports approved CLSI standards</li> </ul>	<p>—/automation SW feature                  automation SW feature/—                  —/—/automation SW feature                  —                  automation SW feature/automation SW feature/—                  automation SW feature                  automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/—                  —/automation SW feature</p>	<p>LIS feature/automation SW feature                  automation SW feature/LIS feature                  automation SW feature/automation SW feature/automation SW feature                  automation SW feature                  automation SW feature/automation SW feature/automation SW feature                  automation SW feature                  automation SW feature                  automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature                  automation SW feature/automation SW feature</p>
<p><b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b></p>	<p>—</p>	<p>Siemens, Cerner, Meditech, SCC Soft, Misys, Data Innovations, OSI, Telepath-iSoft, Netlab, LMX Labzis II, SCL 2000, others/ASTM</p>
<p><b>Transportation systems available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions* (H x W x D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Containers device accommodates/Average throughput in cm per second</li> <li>• Supports automatic rerouting for reflex-repeat-dilutions</li> <li>• Modular HW/Installed options/Device can operate in track and manual mode</li> <li>• Required utilities/Required maintenance</li> <li>• Carrier type/Scalable system</li> </ul>	<p>no                  —                  —                  —                  —                  —</p>	<p>yes                  —/950 x 2,000 x 530 mm/yes                  16, 13 x 100; 16, 13 x 75, others/71.6                  yes                  yes/floor and subfloor mounted/yes                  compressed air, electricity, water/weekly, monthly, quarterly, annually                  single specimen container per carrier/yes</p>
<p><b>Automated centrifugation available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimensions (H x W x D)/Conforms to CLSI Stand. Auto 1-5</li> <li>• Maximum throughput/Containers device accommodates</li> <li>• Can identify tube types for custom programmed rate and spin times per run</li> <li>• More than one centrifuge can be connected to track system</li> <li>• For multi-unit centrifuge, each centrifuge operates independently for rate and time</li> <li>• Maintenance required</li> </ul> <p><b>Automated input/accessioning available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Dedicated lanes for stat samples</li> <li>• Maximum No. of samples that can be loaded/Maintenance required</li> </ul> <p><b>Automated decapping available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> <li>• Removes multiple size tube caps per run/Removes screw type sample caps</li> </ul> <p><b>Automated sorting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Software can sort by</li> </ul> <p><b>Specimen integrity monitor available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> <li>• Containers device accommodates/Maintenance required</li> </ul> <p><b>Automated aliquotting available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput**</li> </ul> <ul style="list-style-type: none"> <li>• Containers device accommodates</li> <li>• Inspects samples for bar code/Detects and reports clots in specimen</li> <li>• Detects and reports quantity not sufficient specimens/Maintenance required</li> </ul>	<p>no                  —                  —                  —                  —                  —                  yes                  —/—/yes/1200                  16, 13 x 100; 16, 13 x 75; 13 x 65 to 16 x 100/—                  up to 600, configuration-dependent/quarterly                  yes                  —/configuration-dependent/yes/1,200                  16, 13 x 100; 16, 13 x 75, multiple/quarterly                  yes/yes                  yes                  —/configuration-dependent/yes/1,200                  16, 13 x 100; 16, 13 x 75, multiple/specimen, method, output                  yes                  —/configuration-dependent/yes/700                  16, 13 x 100; 16, 13 x 75, multiple/quarterly                  yes                  —/configuration-dependent/yes/dependent upon number of aliquots and their volumes                  16, 13 x 100; 16, 13 x 75, multiple                  yes/yes                  yes/quarterly</p>	<p>yes                  —/1,900 x 1,570 x 860 mm/yes                  300/16, 13 x 100; 16, 13 x 75, others                  yes                  yes                  yes                  weekly, monthly, quarterly, annually                  yes                  sample manager/1,900 x 2,040 x 860 mm/yes/325                  16, 13 x 100; 16, 13 x 75, others/yes                  1,000/weekly, monthly, quarterly, annually                  yes                  —/included in centrifuge module/yes/240; independent module/550                  16, 13 x 100; 16, 13 x 75, others/weekly, monthly, quarterly, annually                  yes/yes                  yes                  sample manager/1,900 x 2,040 x 860 mm/yes/325                  16, 13 x 100; 16, 13 x 75, others/specimen, method, output                  onboard each instrument                  integrated on chemistry instrument                  16, 13 x 100; 16, 13 x 75, others/—                  no                  —                  —                  —</p>
<p><b>Instrument (analyzer) interfaces</b></p> <ul style="list-style-type: none"> <li>• Rules-based instrument interface control subsystem</li> <li>• Process control of instrument via control subsystem</li> </ul> <p><b>Physical/hardware (instrument/specimen) interface</b></p> <ul style="list-style-type: none"> <li>• Hematology/Chemistry/Coagulation</li> <li>• Immunoassay/Urinalysis</li> </ul>	<p>no                  no                  —                  —</p>	<p>yes                  yes                  robotic arm interface/point-of-reference sampling/robotic arm interface                  point-of-reference sampling, robotic arm interface/point-of-reference sampling</p>
<p><b>Instruments to which your system or product is interfaced</b></p> <p><b>Other robotic products/components to which system or product is linked</b></p>	<p>—                  —</p>	<p>ADVIA 120/2120i/Autoslide solution, ADVIA 1800/2400 solution, ADVIA Centaur XP solution, CLINITEK Atlas solution, Dimension RxL Max solution, others                  —</p>
<p><b>Automated recapper or sealer available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput*</li> <li>• Recaps-seals multiple size tubes simul./Containers device accommodates</li> <li>• Maintenance required</li> </ul>	<p>recapper                  —/configuration-dependent/yes/1,200                  yes/16, 13 x 100; 16, 13 x 75; 13-16 mm in diameter                  quarterly</p>	<p>no                  —                  —                  —</p>
<p><b>Automated storage and retrieval available</b></p> <ul style="list-style-type: none"> <li>• Model/Dimen. (H x W x D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput*</li> <li>• Containers device accommodates/Connects to the track</li> <li>• Room temperature/Minimum and maximum number of tubes stored per module</li> <li>• Multiple size tubes can be stored in the same module/Maintenance required</li> <li>• Refrigerated storage and retrieval capability</li> </ul> <p><b>Longitudinal upgrade pathway or plan to protect users' investments</b></p>	<p>no                  —                  —                  —                  —                  systems are upgradable</p>	<p>yes                  sample manager/1,900 x 2,040 x 860 mm/yes/325                  16, 13 x 100; 16, 13 x 75, others/yes                  yes/1 and 1,000                  yes/weekly, monthly, quarterly, annually                  no                  flexible and expandable: can contain as few as 2 interfaced components-instruments and can expand to up to 16 interfaced configuration dependent/Siemens Healthcare Diagnostics/24-7                  no/yes</p>
<p><b>Average time to install/Who provides service, support/Hours support is available</b>  <b>On-site biomedical engineer required/User group meets regularly</b></p>	<p>2-3 weeks/Sarstedt or authorized Sarstedt service company/contract dependent                  no/no</p>	<p></p>
<p><b>List price</b></p> <p><b>Individual list prices for components</b></p> <ul style="list-style-type: none"> <li>• Process control SW/Transportation systems/Auto. centrifugation</li> <li>• Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval</li> <li>• Specimen integrity monitor/Automated aliquot</li> <li>• Instrument (analyzer) interfaces/Automated recap</li> </ul>	<p>varies                  —                  —                  —                  —</p>	<p>varies                  —                  —                  —                  —</p>
<p><b>Distinguishing features (supplied by company)</b>                  * For basic building block unit                  ** Average throughput in specimen containers per hour per device</p>	<p>bulk loading module: tubes are dumped into a hopper, eliminating need for pre-racking; modular design enables configuration based on individual requirements; screw-cap recapping; manufacturer of instr. and corresponding consumables</p>	<p>high-throughput lab automation with broad menu, single LIS connection, flexible configurations and ability to connect multiple disciplines with same track system: (chemistry, immunoassay, hematology, coagulation, urine)</p>

Tabulation does not represent an endorsement by the College of American Pathologists.



## Laboratory automation systems and workcells

<b>Part 12 of 14</b>	Siemens Healthcare Diagnostics Tim Keating 511 Benedict Avenue, Tarrytown, NY 10591 302-631-9482 www.usa.siemens.com/diagnostics	Siemens Healthcare Diagnostics Tim Keating 511 Benedict Avenue, Tarrytown, NY 10591 302-631-9482 www.usa.siemens.com/diagnostics
See <a href="http://captodayonline.com/productguides">captodayonline.com/productguides</a> for an interactive version of guide		
<b>Name of system/First year installed/No. of 2011 contracts signed No. of live sites installed in N. America/Europe/Asia-Australia</b>	StreamLab Analytical Workcell/2002/— >160 U.S./>295 worldwide	VersaCell System/2002/80 160/>1,000 worldwide
<b>Automation products that are available</b> • Pre-analytical processor/Total laboratory automation • Automated functions: Accessioning/Track load/Centrifugation/Decapping • Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing • Automated functions: Storage-retrieval/Intelligent sample routing • SW: Dedicated Process Control/Middleware control using LIS/Architecture • Company has dedicated automation support team/Remote system monitoring	yes/yes yes/yes/yes/yes yes/analyzer removes aliquot/no/yes yes/yes yes/yes/open yes/yes	yes/no no/no/no/no yes/no/no/no yes/yes yes/yes/closed yes/yes
<b>Software features/functionality</b> • Patient demographics and insurance data/Rules-based architecture • Supports data retrieval/Internet connectivity • Online real-time help system/QC/Stats and management reports • Evaluates validity and releasability of results from automated analyzers • Specimen tracking/Priority processing/Random-access spec. movement • Supports accession number redundancy (duplicate specimen ID) • Supports specimen carrier and level identification • Unique bar-code number per container required • Specimen routing/Multistop routing (one tube to multiple workstations) • Specimen scheduling/Instrument scheduling • Routes test to workstation/Automatic reflex, repeat, dilutions • Supports multiple HW configuration/Supports other proprietary transport. HW • Sample storage and retrieval SW/Supports approved CLSI standards	automation SW and LIS feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature/automation SW feature automation SW feature automation SW feature/automation SW feature/automation SW feature automation SW feature automation SW and LIS feature automation SW and LIS feature automation SW feature/automation SW feature automation SW and LIS feature automation SW and LIS feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature	LIS feature/automation SW feature automation/SW feature/— —/LIS feature/automation SW feature LIS feature automation SW feature/automation SW feature/automation SW feature automation SW feature automation SW and LIS feature automation SW and LIS feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/—
<b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b>	Cerner, Meditech, SCC, Misys, CHCS, LabGem, Swiss Lab, Medicom, IZASA, Confidentia, others/DBASTM, Dimension Protocol, HL7, ASTM	—/ASTM
<b>Transportation systems available</b> • Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Containers device accommodates/Average throughput in cm per second • Supports automatic rerouting for reflex-repeat-dilutions • Modular HW/Installed options/Device can operate in track and manual mode • Required utilities/Required maintenance • Carrier type/Scalable system	yes StreamLab/60 × 70 × 35 inches/yes 16, 13 × 100; 16, 13 × 75/300 tubes per hour yes yes/floor mounted/yes compressed air, electricity/weekly single specimen container per carrier/yes	yes VersaCell System/70 × 51 × 41 inches/— 16, 13 × 100; 16, 13 × 75/— yes yes/floor mounted/yes electricity/annually single specimen container per carrier/yes
<b>Automated centrifugation available</b> • Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Maximum throughput/Containers device accommodates • Can identify tube types for custom programmed rate and spin times per run • More than one centrifuge can be connected to track system • For multi-unit centrifuge, each centrifuge operates independently for rate and time • Maintenance required	yes StreamLab/31 × 23 × 29 inches/yes up to 400 per hour/16, 13 × 100; 16, 13 × 75, handles various sizes simultan. yes no — weekly, monthly	no — — — — —
<b>Automated input/accessioning available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Dedicated lanes for stat samples • Maximum No. of samples that can be loaded/Maintenance required	yes StreamLab/60 × 70 × 35 inches/yes/300 tubes 16, 13 × 100; 16, 13 × 75/yes up to 600/daily, monthly	yes VersaCell System/70 × 51 × 41 inches/—/200 16, 13 × 100; 16, 13 × 75/yes 200/annually
<b>Automated decapping available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required • Removes multiple size tube caps per run/Removes screw type sample caps	yes StreamLab/integrated with input-output track/yes/300 16, 13 × 100; 16, 13 × 75/daily, monthly yes/yes	no — —
<b>Automated sorting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Software can sort by	yes StreamLab/integrated with input-output track/yes/300 16, 13 × 100; 16, 13 × 75/specimen, method, output	yes VersaCell System/70 × 51 × 41 inches/no/200 16, 13 × 100; 16, 13 × 75/—
<b>Specimen integrity monitor available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required	yes StreamLab/integrated with analyzer/yes/300 16, 13 × 100; 16, 13 × 75/—	no — —
<b>Automated aliquotting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates • Inspects samples for bar code/Detects and reports clots in specimen • Detects and reports quantity not sufficient specimens/Maintenance required	yes StreamLab/integrated with sample transfer module/yes/300 16, 13 × 100; 16, 13 × 75 yes/yes yes/daily	no — — — —
<b>Instrument (analyzer) interfaces</b> • Rules-based instrument interface control subsystem • Process control of instrument via control subsystem Physical/hardware (instrument/specimen) interface • Hematology/Chemistry/Coagulation • Immunoassay/Urinalysis	yes yes no/pt-of-ref sampling, rob. arm interface/pt-of-ref sampling, rob. arm interf. point-of-reference sampling, robotic arm interface/no	yes yes no/point-of-reference sampling/no point-of-reference sampling/no
<b>Instruments to which your system or product is interfaced</b>	Dimension RxL Max, Dimension Vista 1500/500, Immulite 2000 and 2500; Sysmex CA 7000; Dimension EXL with LM, Advia Centaur	Advia 1800, Immulite Immunoassay, Advia Centaur, Dimension EXL with LM, Dimension EXL 200, Dimension RxL MAX
<b>Other robotic products/components to which system or product is linked</b>	—	StreamLab analytical workcell and Advia automation workcells
<b>Automated recapper or sealer available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Recaps-seals multiple size tubes simult./Containers device accommodates • Maintenance required	yes StreamLab/40 × 36 × 17 inches/yes/300 yes/13 × 100; 13 × 75; 16 × 100; 16 × 75 daily, monthly	no — — —
<b>Automated storage and retrieval available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Connects to the track • Room temperature/Minimum and maximum number of tubes stored per module • Multiple size tubes can be stored in the same module/Maintenance required • Refrigerated storage and retrieval capability Longitudinal upgrade pathway or plan to protect users' investments	yes StreamLab SW and input-output module/60 × 70 × 35 inches/yes/300 13 × 100; 13 × 75; 16 × 100; 16 × 75 (47,952 storage capacity)/no yes/up to 576 yes/— yes StreamLab systems are scalable with open configurations	no — — — — continue connectivity development and software enhancements
<b>Average time to install/Who provides service, support/Hours support is available On-site biomedical engineer required/User group meets regularly</b>	five days/Siemens/24–7 no/yes	two days/Siemens Healthcare Diagnostics/24–7 no/no
<b>List price</b> Individual list prices for components • Process control SW/Transportation systems/Auto. centrifugation • Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval • Specimen integrity monitor/Automated aliquot • Instrument (analyzer) interfaces/Automated recap	— — — — —	— — — — —
<b>Distinguishing features (supplied by company)</b> * For basic building block unit ** Average throughput in specimen containers per hour per device Note: a dash in lieu of an answer means company did not answer question or question is not applicable	integrated automation solution with open architecture allows custom configuration and reconfiguration by incorporating a 90-degree track turn, which helps maintain a small footprint	breadth of menu with flexibility of connectivity; throughput, pre and postanalytical sample management

## Laboratory automation systems and workcells

<b>Part 13 of 14</b>	<b>Sysmex America</b> Nilam Patel pateln@sysmex.com 1 Nelson C. White Parkway, Mundelein, IL 60060 800-379-7639 ext. 4309 www.sysmex.com/automation	<b>Sysmex America</b> Krista Curcio curciok@sysmex.com 1 Nelson C. White Parkway, Mundelein, IL 60060 800-379-7639 ext. 4613 www.sysmex.com/automation
See captodayonline.com/productguides for an interactive version of guide		
<b>Name of system/First year installed/No. of 2011 contracts signed</b> <b>No. of live sites installed in N. America/Europe/Asia-Australia</b>	<b>HST-N/1991/50+</b> 350/1,600+ (Europe, Asia, Latin America, Canada, & Australia)	<b>XE-Alpha N/1991/30</b> 250/650+ (Europe, Asia, Latin America, Canada, Australia)
<b>Automation products that are available</b> • Pre-analytical processor/Total laboratory automation • Automated functions: Accessioning/Track load/Centrifugation/Decapping • Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing • Automated functions: Storage-retrieval/Intelligent sample routing • SW: Dedicated Process Control/Middleware control using LIS/Architecture • Company has dedicated automation support team/Remote system monitoring	no/no yes/no/no/no yes/no/—/no no/yes yes/yes/closed yes/yes	—/— yes/—/no/no yes/no/—/no no/— yes/yes/closed yes/yes
<b>Software features/functionality</b> • Patient demographics and insurance data/Rules-based architecture • Supports data retrieval/Internet connectivity • Online real-time help system/QC/Stats and management reports • Evaluates validity and releasability of results from automated analyzers • Specimen tracking/Priority processing/Random-access spec. movement • Supports accession number redundancy (duplicate specimen ID) • Supports specimen carrier and level identification • Unique bar-code number per container required • Specimen routing/Multistop routing (one tube to multiple workstations) • Specimen scheduling/Instrument scheduling • Routes test to workstation/Automatic reflex, repeat, dilutions • Supports multiple HW configuration/Supports other proprietary transport. HW • Sample storage and retrieval SW/Supports approved CLSI standards	automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature /automation SW feature /LIS feature automation SW feature automation SW feature/automation SW feature/yes automation SW feature automation SW feature automation SW feature/automation SW feature —/— automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature/automation SW feature	—/automation SW feature automation SW feature/LIS feature automation SW feature /automation SW feature /LIS feature automation SW feature automation SW feature/automation SW feature/— automation SW feature automation SW feature automation SW feature automation SW feature/automation SW feature —/— automation SW feature/automation SW feature —/automation SW feature —/—
<b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b>	Cerner (Classic and Millennium), Misys, SCC, Meditech, GE/HL7 and ASTM	Cerner (Classic and Millennium), Misys, SCC, Meditech, GE/HL7 and ASTM
<b>Transportation systems available</b> • Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Containers device accommodates/Average throughput in cm per second • Supports automatic rerouting for reflex-repeat-dilutions • Modular HW/Installed options/Device can operate in track and manual mode • Required utilities/Required maintenance • Carrier type/Scalable system	yes HST-N/configuration-dependent/yes 16 × 75; 13 × 75/minutes throughput 150/hour; max as high as lab needs/hour yes yes/floor mounted/yes — rack/yes	yes Alpha N/2 × 7.3 × 3.4 feet 16 × 75; 13 × 75/based on number of analyzers no yes/—/yes — rack/no
<b>Automated centrifugation available</b> • Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Maximum throughput/Containers device accommodates • Can identify tube types for custom programmed rate and spin times per run • More than one centrifuge can be connected to track system • For multi-unit centrifuge, each centrifuge operates independently for rate and time • Maintenance required	no — — — — —	no — — — — —
<b>Automated input/accessioning available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Dedicated lanes for stat samples • Maximum No. of samples that can be loaded/Maintenance required	yes — — 200 samples per input module/—	yes — —/no 100 samples per input module/—
<b>Automated decapping available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required • Removes multiple size tube caps per run/Removes screw type sample caps	no — — —	no — — —
<b>Automated sorting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Software can sort by	yes PVT TS series: low-mid volume ~5 × 3 feet; high volume ~6 × 5 feet 13 × 75/specimen, method, output yes (located within the analyzers)	no —/—/yes/— —/— yes (located within the analyzers)
<b>Specimen integrity monitor available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required	— — —	— — —
<b>Automated aliquotting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates • Inspects samples for bar code/Detects and reports clots in specimen • Detects and reports quantity not sufficient specimens/Maintenance required	no — — — —	no — — — —
<b>Instrument (analyzer) interfaces</b> • Rules-based instrument interface control subsystem • Process control of instrument via control subsystem Physical/hardware (instrument/specimen) interface • Hematology/Chemistry/Coagulation • Immunoassay/Urinalysis	yes yes — point-of-reference sampling/—/— —	yes yes — — —
<b>Instruments to which your system or product is interfaced</b>	Bio-Rad Variant II Turbo Link A1C analyzer	—
<b>Other robotic products/components to which system or product is linked</b>	Thermo automation, Lab Interlink/Labotix, IDS	—
<b>Automated recapper or sealer available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput* • Recaps-seals multiple size tubes simult./Containers device accommodates • Maintenance required	no — — —	no — — —
<b>Automated storage and retrieval available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput* • Containers device accommodates/Connects to the track • Room temperature/Minimum and maximum number of tubes stored per module • Multiple size tubes can be stored in the same module/Maintenance required • Refrigerated storage and retrieval capability Longitudinal upgrade pathway or plan to protect users' investments	no — — — — — —	no — — — — — —
<b>Average time to install/Who provides service, support/Hours support is available</b> On-site biomedical engineer required/User group meets regularly	<3 days/Sysmex/24-7 no/no	1 day/Sysmex/24-7 no/no
<b>List price</b> Individual list prices for components • Process control SW/Transportation systems/Auto. centrifugation • Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval • Specimen integrity monitor/Automated aliquot • Instrument (analyzer) interfaces/Automated recap	varies — — — —	varies — — — —
<b>Distinguishing features (supplied by company)</b> * For basic building block unit ** Average throughput in specimen containers per hour per device Note: a dash in lieu of an answer means company did not answer question or question is not applicable	scalable, flexible, and reliable automation and instrument systems; fast installation (<3 days); scalable multi-site, multi-system middleware solutions that are developed, tested, and supported by Sysmex	scalable and flexible configurations; one-day installation; scalable middleware solutions are developed and supported by Sysmex

## Laboratory automation systems and workcells

<b>Part 14 of 14</b>	<b>Yaskawa America Inc., Motoman Robotics Division</b> Craig Rubenstein craig.rubenstein@motoman.com 100 Automation Way, Miamisburg, OH 45342 949-263-2648 www.motoman.com/labauto/	<b>Yaskawa America Inc., Motoman Robotics Division</b> Craig Rubenstein craig.rubenstein@motoman.com 100 Automation Way, Miamisburg, OH 45342 949-263-2648 www.motoman.com/labauto/
See captodayonline.com/productguides for an interactive version of guide		
<b>Name of system/First year installed/No. of 2011 contracts signed</b> <b>No. of live sites installed in N. America/Europe/Asia-Australia</b>	<b>Autosorter II/2006/4</b> 22/—/—	<b>Autosorter III/2008/4</b> 22/—/—
<b>Automation products that are available</b> • Pre-analytical processor/Total laboratory automation • Automated functions: Accessioning/Track load/Centrifugation/Decapping • Automated functions: Rack specific sort/Aliquot/Tube relabeling/Resealing • Automated functions: Storage-retrieval/Intelligent sample routing • SW: Dedicated Process Control/Middleware control using LIS/Architecture • Company has dedicated automation support team/Remote system monitoring	yes/no yes/yes/yes/yes yes/yes/no/yes (recapping) no/yes yes/yes/open yes/yes	yes/no yes/yes/yes/yes yes/yes/no/yes (recapping) no/yes yes/yes/open yes/yes
<b>Software features/functionality</b> • Patient demographics and insurance data/Rules-based architecture • Supports data retrieval/Internet connectivity • Online real-time help system/QC/Stats and management reports • Evaluates validity and releasability of results from automated analyzers • Specimen tracking/Priority processing/Random-access spec. movement • Supports accession number redundancy (duplicate specimen ID) • Supports specimen carrier and level identification • Unique bar-code number per container required • Specimen routing/Multistop routing (one tube to multiple workstations) • Specimen scheduling/Instrument scheduling • Routes test to workstation/Automatic reflex, repeat, dilutions • Supports multiple HW configuration/Supports other proprietary transport. HW • Sample storage and retrieval SW/Supports approved CLSI standards	—/automation SW feature automation SW feature/automation SW feature automation SW feature/ automation SW feature/automation SW feature — automation SW feature/ automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature automation SW feature automation SW feature/ automation SW feature —/— automation SW feature/— automation SW feature/automation SW feature automation SW feature/automation SW feature	—/automation SW feature automation SW feature/automation SW feature automation SW feature/ automation SW feature/automation SW feature — automation SW feature/ automation SW feature/automation SW feature automation SW feature/automation SW feature automation SW feature automation SW feature automation SW feature/ automation SW feature —/— automation SW feature/— automation SW feature/automation SW feature automation SW feature/automation SW feature
<b>LIS(s) and versions interfaced and live w/LAS/How LIS(s) are interfaced with your LAS</b>	Cerner, Triple G, Surround/ODBC, HL7	Cerner, Triple G, Surround/ODBC, HL7
<b>Transportation systems available</b> • Model/Dimensions* (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Containers device accommodates/Average throughput in cm per second • Supports automatic rerouting for reflex-repeat-dilutions • Modular HW/Installed options/Device can operate in track and manual mode • Required utilities/Required maintenance • Carrier type/Scalable system	yes —/configuration-dependent/yes 16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height/50 yes no/floor mounted/yes compressed air, electricity/daily, monthly, annually single and multiple (30) specimen container per carrier/yes	yes —/configuration-dependent/yes 16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height/50 no no/floor mounted/yes electricity/daily, monthly, annually single specimen container per carrier/yes
<b>Automated centrifugation available</b> • Model/Dimensions (H × W × D)/Conforms to CLSI Stand. Auto 1-5 • Maximum throughput/Containers device accommodates • Can identify tube types for custom programmed rate and spin times per run • More than one centrifuge can be connected to track system • For multi-unit centrifuge, each centrifuge operates independently for rate and time • Maintenance required	no — — — — —	yes Hettich Rotanta/81 × 87 × 42 inches, 9–16 mm diameter, 75–100 mm height/yes 300+/16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height no no — daily, monthly, annually
<b>Automated input/accessioning available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Dedicated lanes for stat samples • Maximum No. of samples that can be loaded/Maintenance required	yes AutoSorter II/6 × 5 × 5 feet/yes/1,000 16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height/yes 1,000/daily, monthly, annually	yes AutoSorter III/81 × 87 × 42 inches (enclosed within ASIII footprint)/yes/800 16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height/yes 300/daily, monthly, annually
<b>Automated decapping available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required	yes —/fits within footprint of AutoSorter II/yes/1,000 16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height/daily, monthly, annually	yes AutoSorter III/81 × 87 × 42 inches (enclosed within ASIII footprint)/yes/800 16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height/daily, monthly, annually
• Removes multiple size tube caps per run/Removes screw type sample caps	yes/yes	yes/yes
<b>Automated sorting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Software can sort by	yes AutoSorter II/6 × 5 × 5 feet/yes/1,000 16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height/specimen, method, output	yes AutoSorter III/81 × 87 × 42 inches/yes/800 16, 13 × 100; 16, 13 × 75, 9–16 mm diameter, 75–100 mm height/specimen, method, output
<b>Specimen integrity monitor available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates/Maintenance required	— — —	— — —
<b>Automated aliquotting available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput** • Containers device accommodates • Inspects samples for bar code/Detects and reports clots in specimen • Detects and reports quantity not sufficient specimens/Maintenance required	yes Aloka APS/68 × 101 × 43 inches/yes/500 16, 13 × 100; 16, 13 × 75 yes/yes yes/daily, monthly, annually	planned Aloka module/to be determined/yes/100–200 16, 13 × 100; 16, 13 × 75 yes/yes yes/daily, monthly, annually
<b>Instrument (analyzer) interfaces</b> • Rules-based instrument interface control subsystem • Process control of instrument via control subsystem	no no	no no
<b>Physical/hardware (instrument/specimen) interface</b> • Hematology/Chemistry/Coagulation • Immunoassay/Urinalysis	Sysmex HST —	Sysmex HST —
<b>Instruments to which your system or product is interfaced</b> <b>Other robotic products/components to which system or product is linked</b>	Sysmex hematology automation MDS (now Innotek) single-specimen carrier transportation system	— ILAS, MDS (now Innotek) single-specimen carrier transportation system
<b>Automated recapper or sealer available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput* • Recaps-seals multiple size tubes simult./Containers device accommodates • Maintenance required	yes (recapper) AutoSorter II/6 × 5 × 5 feet/yes/>1,800 yes/16, 13 × 100; 16, 13 × 75 daily, monthly, annually	planned AutoSorter III/to be determined/yes/800 yes/16, 13 × 100; 16, 13 × 75 daily, monthly, annually
<b>Automated storage and retrieval available</b> • Model/Dimen. (H × W × D)/Conforms to CLSI Stand. Auto 1-5/Avg. throughput* • Containers device accommodates/Connects to the track • Room temperature/Minimum and maximum number of tubes stored per module • Multiple size tubes can be stored in the same module/Maintenance required • Refrigerated storage and retrieval capability	yes — — — — —	yes — — — — —
<b>Longitudinal upgrade pathway or plan to protect users' investments</b>	flexible, open design permits change of tubes/racks as instrumentation changes; connectivity and functionality upgrades	flexible, open design permits change of tubes/racks as instrumentation changes; connectivity and functionality upgrades
<b>Average time to install/Who provides service, support/Hours support is available</b> <b>On-site biomedical engineer required/User group meets regularly</b>	<1–2 weeks, more for complex systems/Motoman/24–7 no/no	<1 week/Motoman/24–7 no/no
<b>List price</b> <b>Individual list prices for components</b> • Process control SW/Transportation systems/Auto. centrifugation • Auto. input, accession/Auto. decap/Auto. sort/Auto. storage and retrieval • Specimen integrity monitor/Automated aliquot • Instrument (analyzer) interfaces/Automated recap	\$250,000 included/configuration-dependent/— included/configuration-dependent/included/— —/configuration-dependent —/configuration-dependent	\$195,000 included/configuration dependent/\$39,500 included/included/included/— —/to be determined —/to be determined
<b>Distinguishing features (supplied by company)</b> * For basic building block unit ** Average throughput in specimen containers per hour per device <i>Note: a dash in lieu of an answer means company did not answer question or question is not applicable</i>	customization-friendly; designed and built in the U.S.; independent of IVD instrument manufacturers; free-standing, high-throughput instruments or integrated lines	customization-friendly; designed and built in the U.S.; independent of IVD instrument manufacturers; free-standing, small footprint, modular automation