

Analytical and precision balances, 9/14

At the AACC Show 2014

September 2014—Ohaus will launch in September its Adventurer line of analytical and precision balances for standard laboratory weighing. Adventurer has nine application modes, including weighing, percent weighing, counting, check weighing, dynamic weighing, totalization, formulation, density determination, and display hold.

The balance can be set up using Adventurer's SmarText 2.0 color touchscreen. Equipped with a virtual QWERTY keypad, the screen offers direct access to operating the balances, inputting user profiles, such as sample, user, or project name, as well as good laboratory practices and other application data. Six mechanical keys provide tactile feedback and one-touch operation of repetitive processes such as tare, zero, calibrate, and print.

The Adventurer is equipped with draft-shield doors constructed of two glass panels, standard connectivity options, including front and back USB ports, and a fully automatic internal calibration system, AutoCal, on select models.

[Ohaus](#), 973-377-9000