

Anti-PIN5 cocktail, 9/13:103

Diagnostic BioSystems' anti-PIN5 cocktail helps distinguish prostate cancer from benign conditions that mimic cancer. The cocktail works through the combined detection of p504S, p40, and high molecular-weight cytokeratin markers. The markers are particularly relevant in helping to diagnose prostatic intraepithelial neoplasia. The markers also offer utility in difficult cases where tissue may be limited.

The anti-PIN5 cocktail is available in ready-to-use format.

Diagnostic BioSystems, 888-896-3350