

# D-dimer kit for Siemens Sysmex CS-2500 analyzers, 9/17

**September 2017**—LGC Maine Standards' Validate D-Dimer linearity and calibration verification test kit has received FDA 510(k) clearance for use on Siemens Sysmex CS-2500 System analyzers. Each kit is prepared using the CLSI EP06-A equal delta sample preparation and offers five distinct concentrations covering the reportable range for Siemens Sysmex CS-2500 System analyzers. Kits are liquid and ready to use.

[LGC Maine Standards](#), 207-892-1300