

[Eppendorf adjustable spacing pipettes](#)

written by CAP TODAY

December 17, 2020

December 2020—Eppendorf released its Move It Adjustable Spacing Pipettes. Move It was designed for the synchronous pipetting of a series of samples between different vessel formats, such as between tubes and plates, without tubing connections between cone and piston-cylinder system. The format of microplates can be adjusted directly with the format limiters. For other sample tubes or tubes in holders, the cone spacing is selected by using the rotary knob.

Move It is available with four, six, eight, and 12 channels, mechanical or electronic.

[Eppendorf](#), 800-645-3050



©2026 CAP TODAY, all rights reserved.