

Hematology control, 8/13:89

Streck's Retic-Chex Linearity, designed to establish the patient reticulocyte reportable range and linearity, is now assayed for the Sysmex XN series (XN-10 and XN-20) hematology instruments.

Values are also assigned for Coulter DxH 800, LH 500, LH 750/LH 755, LH 780, Stks/Maxm/HmX; Abbott Cell-Dyn 3200, 3500/3700, Ruby, Sapphire; and Sysmex XE-2100, XT-2000i/XT-4000i.

The product is available in five-level sets with 3.0-mL plastic cap-pierceable vials. Closed-vial stability is 105 days, and open-vial stability is five days.

Streck, 800-843-0912