

Heparin test kit for Stago STA analyzers, 3/18

March 2018—LGC Maine Standards announced the release of its Validate Heparin linearity and calibration verification kit. The Validate Heparin kit evaluates heparin anti-Xa activity in a human plasma matrix and is prepared using CLSI's EP06-A recommended equal delta method for linearity testing. Each kit offers five distinct concentrations covering the 0.1–2.0 IU/mL reportable range of Stago STA analyzers and is liquid, ready to use.

[LGC Maine Standards](#), 207-892-1300