

HRP Green, 8/13:88

Bio SB now offers HRP Green, a vibrant chromogen to help clinical and research laboratories detect IHC markers on paraffin-embedded tissues.

HRP Green comes in a two-part buffer and chromogen format that is compatible with all Bio SB polymer and biotin HRP-based detection systems. HRP Green requires a five-minute incubation period, can be counterstained with Nuclear Fast Red or hematoxylin, and is xylene (PermaMounter) compatible.

To view promotional brochures that highlight HRP Green as well as detection systems and chromogens that Bio SB offers, visit www.biosb.com/product-brochures.