

[Illumina receives EUA for sequencing-based COVID-19 dx](#)

written by CAP TODAY
August 19, 2020

August 2020—Illumina announced that the FDA issued an emergency use authorization for the Illumina COVIDSeq test, a high-throughput, sequencing-based, in vitro diagnostic workflow enabling the detection of SARS-CoV-2.

COVIDSeq uses upper respiratory specimens, including a nasopharyngeal or oropharyngeal swab, and delivers sample receipt to result in 24 hours using the NovaSeq 6000 sequencing system. The differentiated diagnostic design includes 98 amplicons that target the full SARS-CoV-2 genome.

The workflow accommodates up to 3,072 samples per NovaSeq run, leveraging the S4 flow cell, and includes steps for viral RNA extraction, RNA-to-cDNA conversion, PCR, library preparation, sequencing, and report generation. Components include the NovaSeq 6000 coupled with Illumina Tagmentation library preparation technology and the Dragen COVIDSeq Test Pipeline for rapid reporting.

[Illumina](#), 858-202-4566



©2026 CAP TODAY, all rights reserved.