

Linearity FD D-Dimer PathFast Immunoanalyzer, 3/14

[printfriendly]

March 2014—Audit MicroControls added the Linearity FD D-Dimer PathFast Immunoanalyzer, five levels, to its line of calibration verification/linearity products. It is intended to simulate human patient serum samples for the purpose of determining linearity, calibration verification, and verification of reportable range for D-dimer.

Packaged as 5 × 1 mL, this lyophilized product has a shelf life of two years and open-vial stability of two days when stored at 2°C–8°C. Audit provides a free online, real-time data-reduction program that offers Levey-Jennings charts and peer group analysis for daily quality control product users, as well as linearity graphs and peer group analysis.

A FaxBack option allows customers to fax in their calibration verification/linearity data results. In about one business day, Audit technical support will fax back data-reduction charts, graphs, statistics, and peer group data.

[Audit MicroControls](#), 866-252-8348