

Streck adds instrument to linearity control

February 2019—Streck Retic-Chex Linearity for BC is now assayed for the Beckman Coulter UniCel DxH 900 analyzer. Retic-Chex Linearity for BC is designed to establish the reportable range and linearity of analyzers that report the reticulocyte parameter.

The product is packaged in five-vial sets and is available in 2.5 mL and 3.0 mL plastic cap-pierceable vials. Retic-Chex Linearity for BC has an open-vial stability of five days and a closed-vial stability of 105 days.

[Streck](#), 800-228-6090