

Validate verification kits, 4/18

April 2018—LGC Maine Standards released two Validate linearity and calibration verification kits. The Validate Procalcitonin kit for Roche Cobas analyzers is in a human-serum matrix. The Validate HbA1c kit is for Tosoh analyzers and is in a human, whole blood matrix. It evaluates percentage of glycohemoglobin A1c. The kits are liquid and ready to use, prepared using the CLSI EP06-A equal delta sample preparation, and offer five distinct concentrations covering the reportable range.

LGC Maine Standards MSDRx data reduction software is available at no charge for real-time data analysis, or a laboratory can send its data to the company, where a technical specialist will complete the data analysis and return a report within five business days. Peer group comparison is also available.

[LGC Maine Standards](#), 207-892-1300