

CMV DNA linearity panel, 4/14

April 2014—SeraCare Life Sciences is offering a nine-member panel made from serial dilutions of a cultured virus with established reactivity for CMV DNA called the CMV DNA AccuSpan Linearity Panel.



The panel consists of seven positive members representing serial log dilutions of cultured virus in CMV DNA negative diluent, one negative member prepared from the diluent, and one member of diluent to perform additional dilutions as desired.

The results of the CMV test were calibrated to the WHO International CMV Standard for Nucleic Acid Amplification Techniques and show a nearly 1:1 linearity.

[SeraCare Life Sciences](http://www.seracare.com), 800-676-1881