

## Double-staining kit , 11/17

**November 2017**—Vector Laboratories introduced the VectaFluor Duet kit, offering ready-to-use immunofluorescence double labeling. This kit is ideal for use on a range of frozen and formalin-fixed, paraffin-embedded tissue samples, combining primary and secondary antibodies to streamline workflow. VectaFluor Duet includes a secondary reagent to detect mouse and rat primary antibodies with green and red DyLight fluorescent dyes and a normal blocking serum. With the kit, one protocol replaces several and reduces preparation time, in some cases by as much as a day or more.

Vector Labs also offers a recently published immunohistochemistry guide, which provides a compilation of helpful tools for understanding and implementing IHC. This 48-page guide focuses on all aspects of IHC and its implementation, including history of the technique's development and tips on choosing the right reagents for the application. The "Immunohistochemistry Guide" is available in print and online and can be ordered at no charge or downloaded from the company's website.

**[Vector Laboratories](#)**, 800-227-6666