

Inventory management software, 7/17

October 2017—InvMan, an inventory management software solution from Cove Laboratory Software, automates inventory management by providing direct scanning of vendor reagent barcode labels. Using a portable scanner, the vendor 1-D or 2-D barcode is scanned and uploaded to the reagent database in InvMan. Reagents can be tracked by quantity, lot number, and expiration date. Purchase requisitions can be automatically generated based on reagent need. The equipment inventory database tracks the life cycle and all equipment service visits. When logging in to InvMan, an alert prompts users to expiring lot numbers, reagents to be ordered, and instruments requiring a service visit. InvMan supports the GS1 and HIBC barcode format and installs on a Windows network, and client PCs run the application.

[Cove Laboratory Software](#), 707-938-3075