

## Melanoma cocktail antibody, 2/16

**ScyTek Laboratories announced the availability** of Pan Melanoma Cocktail antibody. This cocktail includes Tyrosinase, MART-1, and HMB45 clones. Tyrosinase and MART-1 are known to be co-expressed in the majority of melanomas. Studies thus far have not shown cases where HMB45 was positive but negative for both MART-1 and Tyrosinase. Therefore the combination of HMB45, MART-1, and Tyrosinase make this cocktail more sensitive and a valuable diagnostic tool for metastatic melanoma.

[ScyTek Laboratories](#), 800-729-8350