

SciGene adds to its FISH product line, 11/13

SciGene's new sodium thiocyanate pretreatment reagent, a convenient 1-M liquid formulation of NaSCN for pretreating tissue samples prior to application of nucleic acid probes for cytogenetic assays, joins an expanding line of FISH slide-processing products from SciGene.

The product is available in a 1-L bottle for Coplin jars or staining dishes and a 4-L container with flow control spout for filling baths on the Little Dipper Processor and VP 2000 instruments. No dilution is required, and the reagent can be stored at room temperature.

[SciGene](#), 800-342-2119