

## Stain QC slides, 2/15

**February 2014—Alpha-Tec Systems' Stain Quality Control Slides** are designed to ensure the quality, accuracy, and reliability of stain reagents and staining procedures. The slides provide a known positive and negative control that the laboratory can use to limit the batch-to-batch differences in staining procedures and standardize diagnostic criteria.

Options include acid-fast bacilli, Nocardia, parasitology, and Gram stain quality control slides in three styles—standard slides with only positive and negative controls or slides that include three wells for patient specimens. Slides are individually sealed in a foil pouch, and each package comes with color-coded log sheets to meet accreditation requirements for documentation.

[Alpha-Tec Systems](#), 360-260-2779