

Thermal profiling system, 4/13:88

SciGene Corp.'s Hyb-Ice thermal profiling system performs everyday temperature profiles. The system consists of a programmable heating/cooling base and removable block that heat and cool samples from 4° to 99°C. It accommodates a variety of quick-change blocks to fit nearly any tube or plate.

The Hyb-Ice system provides excellent temperature regulation, uniformity, and ramp rates and performs PCR in standard- or large-volume tubes. A noncontact heated lid controls condensation.

SciGene Corp., 800-342-2119